

Physical Fitness Test

T.K.K.5401 GRIP-D [DIGITAL GRIP DYNAMOMETER]

Grip-D has been used for over 35 years to measure hand grip strength. After 4 measurements which are carried on left and right hand alternately, it can automatically calculate and display the average of the maximum measured value of each hand.

Specifications:

Measuring range:	5.0 to 100 kg
Minimum measurement unit:	0.1 kg
Display:	LCD in 4 digits
Auto power-off:	1 min after operation
Power source:	Two AAA size batteries
Continuous operating time:	Approx. 100 hr
Grip width range:	40 to 70 mm
Dimension:	154 x 235 x 62 mm
Weight:	0.68 kg



T.K.K.5001 GRIP-A [ANALOG GRIP DYNAMOMETER]

Grip-A has been used for over 50 years to measure hand grip strength. It is an analog indication type grip dynamometer with reinforced plastic cover on the front. It has T.K.K.5001 (General type) and T.K.K.5001b (Schoolchild type).

Specifications:

Measuring range:	0 to 100 kg	(General type)
	0 to 50 kg	(Schoolchild type)
Minimum measurement unit:	0.5 kg	
Dimension:	154 x 235 x 59 mm	
Weight:	0.62 kg	



T.K.K.5825 GRIP-A for Child

Grip width can be adjusted down to 20 mm to measure hand grip strength of children accurately.

Specifications:

Measuring range:	0 to 30 kg
Minimum measurement unit:	0.5 kg
Grip width range:	20 to 50 mm
Dimension:	154 x 235 x 59 mm
Weight:	0.62 kg

