



INFLUENZA A/B ANTIGEN TEST WITH CONTROL SWABS

The ulti med Influenza A/B Antigen test is a rapid chromatographic immunoassay for the qualitative detection of Influenza A and B antigens in nasopharyngeal swab specimens. It is intended to aid in the rapid differential diagnosis of Influenza A and B viral infections.

Format | PU:
Kit, Single pouched

Info:
FloqSwabs CE0123 (MDD) - for details see page 48

Article No.	Form	Specimen	Delivery contents
014L470	Test cassette	Nasal swab specimen	20 Test cassettes 20 FloqSwabs 20 Extraction tubes 20 Dropper tips 1 Extraction buffer 1 Workstation 1 Influenza A Positive Control Swab 1 Influenza B Positive Control Swab 1 Influenza A/B Negative Control Swab 1 Package insert

