

# NON-TOUCH Infrared Thermometer JPD-FR202



1S

1 Second  
Reading


**JPD-FR202** Infrared Thermometer measures the body temperature based on the infrared energy emitted from the forehead. Users can quickly get measurement results after positioning properly the temperature probe at forehead.


## FEATURES

  
Professional  
Accuracy

  
10s Auto  
Shut Off

  
Hygienic  
Non-Touch

  
Smart  
Fever Alert

  
Silent Mode

Fever warning



Safe & Hygienic

Multiple use


20 memories








## DISPLAY DESCRIPTION



## QUICK GUIDE

1. Part hair
2. Wipe sweat
3. Point at the center (1-3cm)



## PACKING SPECIFICATION

CTN SIZE(CM)	QTY/CTN	GW(KG)	FCL/20
47 x 34 x 28	50	9.5	560

### ACCESSORIES

Storage bag × 1  
User manual × 1

AAA batteries × 2



## SPECIFICATION

Display	HD LCD Display
Measuring range & accuracy	
Forehead mode	35-42.2°C ± 0.2°C / 95-108.0°F ± 0.4°F
Object mode	0-100.0°C ± 1.0°C / 32-212.0°F ± 2.0°F
Memory storage	20 sets of Measurement data
Temperature units	°C/°F adjustable
Automatic shutdown function	10s ± 1s
Battery	AAAx2