

If there's something to discover, you'll be sure to see it.

The new HEINE DELTA 30 Dermatoscope.



30
YEARS
1st DERMATOSCOPE



The new HEINE DELTA 30 Dermatoscope. You could buy it because of its looks. Or because of its optics.

Extremely sharp image with the new achromatic HEINE optical system

Fast change from contact to non-contact examination

Ergonomic handling due to angled product design



Largest true field of view of 30 mm (lens 32 mm) Ø

LEDHQ
LED NOW IN HEINE QUALITY.

Excellent colour rendering due to LEDHQ

Digital documentation with the exclusive HEINE DERM app or the HEINE Cube System*

• **Crystal clear image**, no matter if you prefer a closer or greater working distance

• **Polarisation** eliminates reflections and glare

• **Brightness control**: 3 levels of light intensity

• **Easily cleaned with spray disinfectant**

• **10x magnification**

• **No cross-contamination**: autoclavable contact plate with scale

• **Dioptre compensation**

• **Switch between polarised and non-polarised illumination mode at the touch of a button** = toggle function. For more information, please visit heine.com/en/Toggle

• **Rechargeable battery with Li-ion technology**:

no memory effect; safely recharge the battery regardless of its state of charge

• **Small contact plate** for difficult-to-access pigmented lesions

• **Charger 30 table charging station**

• **3 possibilities of digital documentation**:

- With the exclusive HEINE Derm App via iPhone** for model 6 or later: Bodymap, patient folder, e-mail function and much more

- Wireless with the HEINE Cube System for up to 10 doctors and HEINE dermatoscopes

- Conventional digital images with many popular smartphones

DELTA 30 with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-230.28.305
Charger 30 table charging station for DELTA 30	X-002.99.212

• **Accessories and further articles under heine.com/en/DELTA30**

* Not available in all countries.

** The Apple iPhone is not included in the scope of delivery. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7, 82211 Herrsching, Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02, E-Mail: info@heine.com
www.heine.com



01/20. A-000.00.228 e



We reserve the right to change specification without notice.
© 2020 HEINE Optotechnik GmbH & Co. KG. All rights reserved.