



# ASPEL

## 308 ABPM

v.101

<p><b>Software</b></p>	<ul style="list-style-type: none"> <li>- cooperation with ASPEL 508 ABPM and HoICARD CR-07 software, automatic recorder configuration.</li> </ul>
<p><b>Functionality</b></p>	<ul style="list-style-type: none"> <li>- possibility of manual programming without using a computer and dedicated software</li> <li>- ability to quickly start registration without using a computer and dedicated software</li> <li>- built-in default, recommended schedule for recording Holter measurements broken down into activity and sleep</li> <li>- battery operation - full functionality</li> <li>- stationary blood pressure monitor mode - powered by USB or battery</li> <li>- full presentation of measurement data from the current test: current cuff pressure (numerical, graphic value - bargraph), exam results (systolic, diastolic pressure, average and HR), all possible errors in the form of warnings and codes</li> <li>- Extensive menu: information data (e.g. device name, patient ID, start date and time of the exam, number of saved exam and availability of free memory, total mileage, additional service data), quick start - quick start of registration with current settings, all saved examination results (ordinal exam number, exam time, systolic, diastolic pressure, average and HR), restoring factory settings, configuring and programming the recorder for automatic operation</li> </ul>

# ASPEL

## 308 ABPM

v.101

<b>Functionality</b>	<ul style="list-style-type: none"> <li>- pressure gauge test function (static pressure measurement from two independent measuring modules)</li> <li>- leak test function (automatic verification of correct operation of the device, tightness of valves, sensors, pump and cuff tightness)</li> <li>- energy saving in the multi-hour recording mode (using deep sleep)</li> </ul>
<b>Battery</b>	- up to 168 hours (up to a week time), memory for 1000 examination
<b>Keyboard</b>	- useful keyboard and separate event button enabling asynchronous measurements (e.g. in the case of deterioration of well-being)
<b>Screen</b>	- transparent, graphic OLED (262,144 colors)
<b>Physical parameters</b>	<ul style="list-style-type: none"> <li>- dimensions: 90x36x93 mm (L x W x H)</li> <li>- weight: &lt;190 g without batteries</li> <li>- housing: molded ABS plastic</li> </ul>
<b>Electrical parameters</b>	<ul style="list-style-type: none"> <li>- power supply: four LR6 (AA) alkaline or rechargeable batteries (AA), USB port: 5 VDC</li> <li>- current consumption: &lt;400 mA - measurement, &lt;70 mA - ready for work</li> </ul>
<b>Pressure measurement</b>	<ul style="list-style-type: none"> <li>- measuring range: adult mode: (0 - 300) mmHg; pediatric mode (0 - 150) mmHg (relative to atmospheric pressure)</li> <li>- pressure sensor: high class with high accuracy and reliability with a compensation circuit</li> <li>- temperature influence compensation: (0 - 50) °C</li> <li>- analog-digital processing: 24 bits, sampling 1000 sps</li> <li>- deflation: linear, adaptive</li> </ul>
<b>Performance parameters</b>	<ul style="list-style-type: none"> <li>- number of measurement periods: 6, definition of activity / sleep</li> <li>- stored data: name, surname, ID, device name</li> <li>- deflation: automatic and adjustable (1000 - 5000) umHg / s, setting resolution 1 umHg / s</li> <li>- pediatric / adult mode of operation</li> <li>- pressure setting: initial and maximum pressure</li> </ul>
<b>Security norms</b>	<ul style="list-style-type: none"> <li>- type of protection against electric shock (EN 60601-1): internally powered device, device class II</li> <li>- application part (EN 60601-1): type CF</li> <li>- medical device class: IIa (rule 10)</li> <li>- clinical measuring accuracy: EN ISO 81060-2</li> <li>- operational safety: EN 60601-1, EN 60601-2-30</li> <li>- electromagnetic compatibility: EN 60601-1-2</li> </ul>

# ASPEL

## 308 ABPM

v.101

<b>Standard equipment</b>	<ul style="list-style-type: none"><li>- cuff for adults - medium size, shoulder circumference (25 - 35) cm</li><li>- pressure hose</li><li>- cuff and pressure hose connector</li><li>- recorder bag with belt</li><li>- LR6 (AA) alkaline batteries - 4 pcs.</li><li>- instructions for use</li><li>- warranty card</li></ul>
<b>Optional equipment</b>	<ul style="list-style-type: none"><li>- pediatric cuff - small size, shoulder circumference (18 - 26) cm</li><li>- cuff for adults - large size, shoulder circumference (33 - 47) cm</li><li>- one-time cuff</li><li>- recorder bag, disposable</li></ul>
<b>Certificates</b>	<ul style="list-style-type: none"><li>- CE 0197</li><li>- EN ISO 13485</li><li>- MDD/93/42/EEC</li></ul>