

PRODUCT COMPARISON:

Bionix® Curette Tip Styles

Our Bionix® Currettes offer a variety of tip styles to match different types of wax and physician preference.

WHAT ARE YOUR CURETTE TIP OPTIONS?

With our varying curette tip options, you'll be able to find the right tool for every cerumen situation. See what sets each tip apart and how our varying styles can help you meet specific patient care needs.



SAFE EAR CURETTES™	9555	FlexLoop® 4mm All-purpose curette for a variety of patients	2240
	6333	CeraSpoon® 4mm Designed for soft wax and patients who produce excessive amounts of cerumen	2250
	2999	ControLoop® 4mm Designed with a “Shapleigh” tip for pulling wax straight out of the ear canal	N/A
	3444	AngleLoop® 4mm Designed with an angled tip for hard-to-get cerumen in older patients	2210
	4111	VersaLoop® 3mm Smooth, flexible tip designed for increased control and improved patient experience	2270
	4888	InfantScoop® 2mm Designed for smaller ears or use through an otoscope	2260
	1222	MicroLoop® 3mm Designed with a rigid construction to remove tough, impacted cerumen	2220



CORRESPONDING LIGHTED EAR CURETTE™

With our Lighted Ear Curette™, you can combine illumination and magnification in one handheld tool. These features make the cerumen easy to see—giving your patients a quick procedure they can be confident in.

