

PRODUCT COMPARISON:

Bionix® OtoClear® Ear Irrigation Kits

Bionix® OtoClear® Ear Irrigation allows healthcare professionals to utilize a soft, flared, 3 stream tip that eliminates the mess, improves patient safety and is compatible with a variety of delivery devices to meet the personal preference of every facility.



THE UNIQUE DESIGN OF OTOCLEAR® EAR IRRIGATION TIPS

- Directs water to the ear canal wall instead of directly at the tympanic membrane
- Prevents backsplash
- Prevents over insertion

7200 OtoClear Ear Irrigation Tips (40/box)
7202 OtoClear Ear Irrigation Tips (20/box)



OtoClear® WaterPik® Kit



OtoClear Tabletop WaterPik Kit



OtoClear SprayWash Kit



OtoClear Aquabot® Kit

	OtoClear® WaterPik® Kit	OtoClear Tabletop WaterPik Kit	OtoClear SprayWash Kit	OtoClear Aquabot® Kit
VOLUME	210ml self-contained water reservoir	1035ml water reservoir	500ml bottle	650ml bottle
VOLUME DELIVERY	Greater compared to manual devices	Greater compared to manual devices	Greater compared to manual devices	Greater compared to manual devices
IRRIGATION TYPE	Pulsed	Pulsed	Intermittent lavage	Continuous stream and pulsed
NUMBER OF SETTINGS	3	6	N/A	N/A
ADDITIONAL FEATURES	Portable. Rechargeable battery.	N/A	Portable. No charge needed. Allows for two-handed procedures.	Portable. No charge needed. Allows for two-handed procedures.
PRODUCT OFFERINGS	7280 OtoClear WaterPik Kit 1 Portable WaterPik, 40 OtoClear Tips, 1 Ear Basin 7245 Portable WaterPik Kit 1 Portable WaterPik, 5 OtoClear Tips, 1 Ear Basin	7240 Tabletop WaterPik Kit 1 Tabletop WaterPik, 5 OtoClear Tips, 1 Ear Basin	7290 OtoClear SprayWash Kit 1 SprayWash Bottle, 20 OtoClear Tips, 1 Ear Basin 7295 OtoClear SprayWash Kit 1 SprayWash Bottle, 5 OtoClear Tips, 1 Ear Basin	7270 OtoClear Aquabot Kit 1 Aquabot Bottle, 20 OtoClear Tips, 1 Ear Basin 7275 OtoClear Aquabot Kit 1 Aquabot Bottle, 5 OtoClear Tips, 1 Ear Basin

INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS