

MIMSAL



# Healthcare Lighting

2021

**MIMSAL**

MIMSAL

CATALOGUE HEALTHCARE LIGHTING 2021

© MIMSAL TRADE S.L.

C. Mollet 17, Pol. Ind. Palou Nord  
08401 Granollers (Barcelona) - Spain  
Tel. +34 930 139 860  
mimsal@mimsal.com  
www.mimsal.com



**MIMSAL**

MIMSAL is a specialized lighting solutions supplier for the healthcare sector. We offer a wide range of luminaires designed to satisfy the diversity of the end-user needs and for each activity sector, by offering technical advice and commercial support to the professionals / prescribers of the healthcare.

Our products are engineered for easy mounting, with modern electronic components and light sources for the best energy efficiency and economy.

Our products and solutions are manufactured and certified in accordance with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices.

MIMSAL complies with the quality certification ISO 9001 and ISO 13485.



# X-RAY VIEWERS AF LED

Extraslim X-Ray Viewers series with a 25 millimeter thickness and made of aluminum frame. This format provides a LED light source for greater illuminance. This format is characterized by its low consumption and lower temperature of the screen.

Available in three formats:

- AF 400 LED: Portable version with desktop stand (optional).
- AF 300 LED: Wall version.
- AF 200 d/c LED: Wall version.



**Desktop stand AF 400 LED**  
Cod. AF-A



*Accessory stainless steel desktop stand for AF 400 LED*

## AF 400 LED

Technical Details	(Cod. AF400LED)
Light source	LED 8 W
Illuminance	Lux 9950 cd/m <sup>2</sup>
Colour temperature	6000 °K
Lamp lifetime	50.000 hours
Illuminated surface	295 x 210 mm
Total dimensions	365 x 280 x 25 mm
Power supply	110/240 V - 50/60 Hz
Frame material	Aluminum
Recommended temperature	from 10 to 40 °C

## AF 300 LED

Technical Details	(Cod. AF300LED)
Light source	LED 14 W
Illuminance	Lux 6300 cd/m <sup>2</sup>
Colour temperature	6000 °K
Lamp lifetime	50.000 hours
Illuminated surface	350 x 420 mm
Total dimensions	420 x 490 x 25 mm
Power supply	110/240 V - 50/60 Hz
Frame material	Aluminum
Recommended temperature	from 10 to 40 °C

## AF 200 d/c LED

Technical Details	(Cod. AF200d/cLED)
Light source	LED 32 W
Illuminance	Lux 6700 cd/m <sup>2</sup>
Colour temperature	6000 °K
Lamp lifetime	50.000 hours
Illuminated surface	420 x 750 mm
Total dimensions	490 x 820 x 25 mm
Power supply	110/240 V - 50/60 Hz
Frame material	Aluminum
Recommended temperature	from 10 to 40 °C



Front view and extraslim profile view detail  
AF 400 LED version



# ACCESSORIES

## WALL B SUPPORT

Wall support.  
Separation between axis-wal: 7cm  
Cod. B0101800



## AH TABLE CLAMP

Table fixing.  
Cod. B010190



## RAIL PLUS BRACKET

To adapt to U.C.I rail.  
Cod. 82026



## TROLLEY STAND

H=92 cm Ø 61 cm.  
4 kg Cod. 09655  
8,8 kg Cod. 09651



## EXTENSION ARM + WALL B SUPPORT

It provides 40 cm rotation from the wall.  
Cod. B010150





**MIMSAL**

MIMSAL TRADE S.L.

C. Mollet 17, Pol. Ind. Palou Nord  
08401 Granollers (Barcelona) - Spain  
mimsal@mimsal.com - www.mimsal.com