

DYONCOVIDAG



DYONMED S.A.
www.dyonmedical.com

EXPORT THIS PAGE IN PDF

REGULATORY STATUS

CE-IVD

TIME TO RESULT

15

ASSAY TARGET

Antigen

TEST FORMAT

Lateral flow assay (strip or cassette)

TECHNOLOGY PRINCIPLE

Antibody-based detection

TARGET ANALYTE

N protein
S protein (S1)
S-protein (type unknown)
S-protein RBD
S-protein S2

LEVEL OF AUTOMATION

Manual

LABORATORY/POINT-OF-CARE

Point-of-care

SELF-TESTING/SELF-COLLECTION

Professional use only

VALIDATED SAMPLE TYPES

Nasopharyngeal swab
Oropharyngeal swab

CONNECTIVITY

No

Manufacturer performance data

SENSITIVITY

91%

SPECIFICITY

99%

DYONMED-S.A

HIDE FILTERS

FILTER THE TEST LIST

Search tests



APPLY FILTERS

MANUFACTURER

DYONMED S.A. (4)


REGION

Europe (4)


COUNTRY

Greece (4)

Assay target



Manufacturer region



RESULTS (4)

CORONA-Ag
DYONMED S.A.

DyonCovid19-S
DYONMED S.A.

DyonCovidAg
DYONMED S.A.

FluCov2
DYONMED S.A.