






# New Product



## Reusable Pressure Infusion Bag

### Product Advantage

Photo	Feature	Benefit
See picture above	Reusable use	Reduce the economic cost of hospitals
	Metal pressure indicator mark	Precise control of pressure and flow
	Large inflatable ball, soft, good elasticity	Easier and more efficient inflation
	Unique design, equipped with Robert clip	Secondary pressure keeping to avoid air leakage, more safe and reliable
	Unique design, equipped with Robert clip	Secondary pressure keeping to avoid air leakage, more safe and reliable
	Transparent nylon mesh material	Not easy to cut and tear, can clearly observe the infusion bag and residual volume, convenient for rapid setting and replacement of infusion bag

## Scope of Application

- ★ Used to continuously pressurize the liquid containing heparin to flush the built-in arterial piezometer.
- ★ Arterial transfusion during surgery and emergency situations.
- ★ In cerebrovascular surgery interventional treatment, in order to prevent blood from flowing back to the coaxial catheter leading to thrombosis and falling off, causing vascular embolism in the body, it is necessary to use an infusion pressurized bag for high-pressure perfusion, and continuous saline drip irrigation in the catheter.
- ★ Emergency fast blood transfusion in the field, battlefield, hospital and other occasions.

## Product Parameter

Capacity	MED-LINK REF NO.	Specification	Packing	Price Code	Photo
500ml	Y000R05	Metal pressure indicator, silicone tube, material: nylon, waterproof	24pcs/box	B6	C
1000ml	Y000R10	Metal pressure indicator, silicone tube, material: nylon, waterproof	24pcs/box	B0	B
3000ml	Y000R30	Metal pressure indicator, silicone tube, material: nylon, waterproof	24pcs/box	C5	A

## 【Related Product Recommendations】 Disposable Pressure Infusion Bag



\*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
 Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China  
 Email: sales@med-linket.com | Website: www.med-linket.com