

## MEDICATION CART

### Summary

<b>Series</b>	SMART
<b>Item</b>	Medication cart
<b>Code</b>	CF101
<b>Dimensions</b>	mm. 620x590x1028H
<b>CND</b>	V0899
<b>RDM</b>	2025732/R

### Description

Configuration:  
4+1 drawers, push handles, sliding working top, n. 3 accessories rails, central locking device by key.

#### STRUCTURE:

Electro-galvanised steel sheet structure, 10/10 thickness, bake painted on epoxy base, acrylic, unscratchable stain-resistant polyurethane paint finish. Easy to clean and disinfect. Perimetral bumper.

#### UPPER PART:

Worktop made of auto-extinguishing technopolymer, easily washable, with raised edge on three sides.

#### DRAWERS:

5 epoxy-painted sheet steel fully-extractable on telescopic guides with anti pull-out stops, with integrated handles, blue colored. Central locking mechanism. Standard drawers configuration:

- nr. 4 drawers dim. 130h mm,
- nr. 1 drawers dim. 200h mm.

#### LEFT SIDE:

Ergonomic push handle and 1 accessories bar holder.

#### RIGHT SIDE:

n. 2 accessories bar holder.

#### BASE:

At the base nr. 4 non-marking rubber wheels Ø 125 mm, each with brake.

