

EMERGENCY CART

Summary

Series	SMART
Item	Emergency cart High
Code	CF304
Dimensions	mm. 620x590x1070H
CND	V0899
RDM	2025813/R



Description

Configuration:

4+1 drawers, push handles, sliding working top, n. 3 accessories rails, central locking device by seal, kit 100 seal, CPR massage board, adjustable IV pole, defibrillator tray, waste bin, lateral basket, gloves holder.

STRUCTURE:

Electro-galvanised steel sheet structure, 10/10 thickness, bake painted on epoxy base, acrylic, unscratchable stain-resistant polyurethane paint finish. Easy to clean and disinfect. Perimetral bumper.

UPPER PART:

Worktop made of auto-extinguishing technopolymer, easily washable, with raised edge on three sides.

DRAWERS:

5 epoxy-painted sheet steel fully-extractable on telescopic guides with anti pull-out stops, with integrated handles, red colored. Central locking mechanism. Standard drawers configuration:

- nr. 4 drawers dim. 130h mm,
- nr. 1 drawers dim. 200h mm.

LEFT SIDE:

Ergonomic push handle and 1 accessories bar holder, lateral basket.

RIGHT SIDE:

n. 2 accessories bar holder, oxygen bottle holder, waste bin with lid 10 lt, gloves holder, oxygen bottle holder.

BACK:

CPR massage board, IV pole, defibrillator tray.

BASE:

At the base nr. 4 non-marking rubber wheels Ø 125 mm, each with brake.

