











You've Earned the Right to Wear the Best



3M™ Littmann® Stethoscopes Product Comparison Chart

	Electronic Stethoscopes		Cardiology Stethoscopes			Classic Stethoscopes			Select/Lightweight Stethoscopes	
	Electronic Stethoscope Model 3200	Electronic Stethoscope Model 3100	Master Cardiology™ Stethoscope	Cardiology S.T.C. (Soft Touch Chestpiece) Stethoscope	Cardiology III™ Stethoscope	Master Classic II™ Stethoscope	Classic II S.E. Stethoscope	Classic II Pediatric and Infant Stethoscope	Select Stethoscope	Lightweight II S.E. Stethoscope
										
Performance (1-10 Scale)	10+	10+	10	9	9	8	7	7	6	6
Main Application Areas	Cardiology/ High Performance Bluetooth®	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Blood Pressure, Limited Physical Assessment	Blood Pressure, Limited Physical Assessment
Chestpiece Design and Technology	Chrome	Chrome	Polished Stainless Steel Single sided	Matte Finish Stainless Steel Single sided with Soft Touch Insert	Machined Stainless Steel Double sided Tunable	Plated and Polished Alloy Single sided	Machined Stainless Steel Double sided	Machined Stainless Steel Double sided	Coated Alloy Single sided	Metal/Resin Composite Double sided
Headset Material	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy
Overall Length	27 inch/69 cm	27 inch/69 cm	22 inch/56 cm 27 inch/69 cm	27 inch/69 cm	22 inch/56 cm 27 inch/69 cm	27 inch/69 cm	28 inch/71 cm	28 inch/71 cm	28 inch/71 cm	28 inch/71 cm
Total Weight	185 grams	185 grams	175, 185 grams	180 grams	175, 180 grams	160 grams	135 grams	105, 95 grams	120 grams	118 grams
Diaphragm Details, Diameter	51.0 mm	51.0 mm	Tunable 1.75 inch/4.4 cm	Tunable 1.88 inch/4.8 cm	Two Tunable Diaphragms 1.7 inch/4.3 cm 1.3 inch/3.3 cm	Tunable 1.75 inch/4.4 cm	Tunable 1.75 inch/4.4 cm 1.25 inch/3.2 cm	1.31 inch/3.3 cm 1 inch/2.5 cm 1.06 inch/2.7 cm 0.75 inch/1.9 cm	Tunable 1.75 inch/4.4 cm	Tunable 2.125 inch/5.4 cm
Warranty	2 Years	2 Years	7 Years	6 Years	5 Years	3 Years	3 Years	3 Years	2 Years	2 Years

3M™ Littmann® Ordering Information and Tube Color Guide



	Black	Burgundy	Navy Blue	Hunter Green	Plum	Chocolate	Gray	Pine Green	Cell Blue	Purple	Raspberry	Caribbean Blue	Orange	Seafoam Green	Red	Lilac	Royal Blue	Ocean Blue	Peach	Pearl Pink	Coral
Model 3200 27 in. (69 cm)	3200BK27	3200BU	3200NB																		
Model 3100 27 in. (69 cm)	3100BK27	3100BU	3100NB	3100HG	3100PL																
Master Cardiology™ 22 in. (56 cm) 27 in. (69 cm)	2159																				
	2160 2161 ¹ 2175 ² 2176 ⁵	2163	2164	2165	2167							2178									
Cardiology S.T.C. 27 in. (69 cm)	4471	4472	4473		4475																
Cardiology III™ 22 in. (56 cm) 27 in. (69 cm)	3127																				
	3128 3131BE ¹ 3128BRS ² 3157SM ⁶ 3152RBW ⁶	3129	3130	3134	3135	3137 3137CPR ³	3136		3146		3148	3138			3140						3149
Master Classic II™ 27 in. (69 cm)	2144L 2141 ¹ 2142G ⁴	2146	2147	2632				2634	2633	2143		2630									
Classic II S.E. 28 in. (71 cm)	2201 2218BE ¹ 2201BRS ² 2827SM ⁵	2211	2205	2208		2820CPR ³	2203	2818	2813	2209	2210 2829RBW ⁶	2206 2823 ³	2812	2814				2819	2822	2817	
Classic II Pediatric 28 in. (71 cm)	2113		2123					2131			2122	2119 2153 ³	2155		2113R		2136		2154		
Classic II Infant 28 in. (71 cm)	2114							2132			2125 2157 ⁶	2124	2179		2114R	2126	2156		2158		
Select 28 in. (71 cm)	2290	2293						2305	2301	2294	2296	2291				2303	2298	2306	2310		
Lightweight II S.E. 28 in. (71 cm)	2450	2451							2454			2452		2455		2453					

Teaching stethoscopes, catalog numbers 2138 and 2139 as well as veterinary stethoscope, catalog number 1392 are also available.



3M Health Care
Littmann® Brand Products
 3M Center, Building 275-4E-01
 St. Paul, MN 55144-1000
 U.S.A.
 1 800 228-3957
www.3m.com/healthcare
www.Littmann.com/us

3M Health Care
3M Canada
 Post Office Box 5757
 London, Ontario N6A 4T1
 Canada
 1 800 364-3577

3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology III, Master Classic II and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Bluetooth® is a registered trademark of Bluetooth SIG. Please recycle. Printed in U.S.A. © 3M 2012. All rights reserved. 70-2008-9603-6 (12.12.25) DPI



- ¹ Black Edition – black finish chestpiece and eartubes.
- ² Brass Edition – brass finish chestpiece and eartubes.
- ³ Copper Edition – copper finish chestpiece and eartubes.
- ⁴ Gold Edition – gold finish chestpiece and eartubes.
- ⁵ Smoke Edition – smoke finish chestpiece and eartubes.
- ⁶ Rainbow Edition – rainbow finish chestpiece and sunshine gold eartubes.