

 **MDSCOPE**[®]

Video Otoscope




GOLDEN PIN
DESIGN AWARD



reddot design award
winner 2011

On-the-spot examination



Instant View

A clear and sharp image can be acquired and displayed momentarily on a built-in screen, with press of a button.

Real Time Interaction

The medical image can be captured and shared simultaneously with patients.

A true improved interaction between the health professionals and the patients.

Electronic Medical Record

A medical image is automatically generated to support most of the stand-alone health information system. Simply plug the cable into the video-out socket to retrieve the image.

Quality & safety standards include:

Medical Electronic Equipment EN60601-1-2, EN60601-1
Medical Device Directive 93/42/EEC
FDA & CE registered and listed
ISO 13485 & GMP certified manufacturer

Product specifications:

Body size : 2.5" width x 8" height
2.4" color LCD display
Light weight: 218 grams (0.48lb)
Battery lifespan: 4 hours (AAx2)

Patents and awards:

Invention & Design Patents Worldwide
Golden pin design award 2011 - Taiwan
Red dot design award winner 2011 - Germany
Good design award 2012 - USA

Special features include:

Image freeze function
Simultaneous image displayed on TV
Digital image saved to Mac / PC system

Portable and Affordable



Detachable & Upgradable camera probes are designed for different exams as needed.



Elite Pack Series



MS101 (ELB)

- Include 1 Standard Camera Probe 45mm L
- RCA Video Cable
- Disposable Specula (24pcs)
- AA Alkaline Battery (2pcs)

Deluxe Zipper Pack Series



MS101(DX3)

- Include 4 Camera Probes- 45 / 75 / 150 / 300mm L
- Image Dongle
- Cable Set
- Disposable Specula (24pcs)
- AA Alkaline Battery (2pcs)

Limited Edition



MS101 (D1)

- Include 3 Camera Probes- 45 / 75 / 150 mm L
- Image Dongle
- Cable Set
- Disposable Specula (24pcs)
- AA Alkaline Battery (2pcs)



Veterinary Endoscope Set



MS101(VXB)

- Include 3 Camera Probes- 75mm/150mm/300mm L
- Image Dongle
- Cable Set
- Alcohol Prep Pad (20pcs)
- AA Alkaline Battery (2pcs)



Accessories



MS101-002T

- Standard Camera Probe 45mm L for Ear Examination



MS101-015

- Disposable Specula for MS101-002T



MS101-102

- Pediatric Camera Probe 45mm L for Ear Examination



MS101-019

- Silicon Specula for MS101-102



MS101-012T

- Long Camera Probe 75mm L for Ear & Oral Examination



MS101-012V

- Long VGA Camera Probe 75mm L for Ear & Oral Examination

MS101-013T

- Flexible Camera Probe 150mm L for Oral Examination



MS101-014T

- Flexible Camera Probe 300mm L for Long Distance Examination



MS101-011

- Image Dongle to Download Image to PC or Mac System



MS101-020

- Resting Cradle for MS101



A Revolutionary Invention

www.mdscope.net

"MDScope" is a registered trademark of **Apple BioMedical Inc.**

✉ sls@applebmi.com

☎ +8862.2999 5505

☎ +8862.2999 6605