



## Features

- Quality and Reliability
- True Zero Start
- Accurate Measurements
- Whisper Quiet
- High Compatibility
- Stress and Pulmonary Function Testing

## Whisper Quiet

CardioMed Treadmills are designed and engineered to be exceptionally quiet making them ideal for any environment whether it is an Olympic Training Center or office setting.

## Zero Start

CardioMed treadmills have a true zero start and then gently increases the speed in manageable 0.1 mph / .16 km/h increments for patient comfort, safety and confidence.

## Importance of Accurate Speed

High torque is essential for accurate speed and testing. CardioMed Treadmills provide uniform digitally controlled high torque entirely over the treadmill's speed range and the patient's weight range, up to a maximum user weight of 500 lbs/228 kgs. The higher torque alleviates the heavy foot fall and slipping of the running belt. Overall the patient has a more enjoyable treadmill experience and the user is confident that the treadmill is performing accurately.

**NEW!**



*Optional 15" Touch-Screen  
with CardioMed Controller*