

D-piece™

Bacterial filter for use with the SteriBreath™ mouthpiece.

The D-piece™ has been specifically designed to be used with the SteriBreath™ mouthpiece.

Restricted flow: the D-piece™ will restrict the patient's expiration flow by narrowing the air's pathway through the D-piece™.

The D-piece™ contains a one-way valve to prevent patients drawing air back from the monitor. An integrated infection control filter removes and traps >99% of airborne bacteria and 96% of viruses¹. Bedfont® recommend that the D-piece™ is replaced at least every 28 days or if visibly soiled.

¹Public Health England. An evaluation of filtration efficiencies against bacterial and viral aerosol challengers report No. 17/001. Salisbury: PublicHealth England; 2017.



Specification

Material	Polypropylene
Item: D-piece	(12 Month Supply)
Item: D-piece-6	(6 Month Supply)

www.covita.net | 800.707.5751 

SteriBreath™ mouthpieces

Single-use mouthpiece for use with the D-piece™

Combining ease of use with gold standard infection control to create the perfect product for use with our breath monitors. The SteriBreath™ mouthpieces now bring you the following benefits:

- **Excellent infection control:** each mouthpiece is individually wrapped and can be located very easily onto the D-piece™ without having to touch the mouthpiece at all
- **Ease of use:** being round, the mouthpiece does not need to be 'popped' open
- **Space saving:** SteriBreath™ are very thin, resulting in a box of 250 taking up very little space and costing less in shipping compared to traditional cardboard mouthpieces.



Specification

Dimensions	32x14.5x12cm (box of 250)
Weight	0.38kg (box of 250)
Material	Polypropylene

Pack of 250 SteriBreath™ mouthpieces

Order code: STERIBREATH-MP

coVita
212 Cottage Grove Ave, Suite C
Santa Barbara, CA 93101
Tel: 800.707.5751 Fax: 800.721.2377
Email: service@covita.net
www.covita.net

© Bedfont® Scientific Limited 2017

Issue 2 - May 2017, Part No: MKT336
Bedfont® Scientific Limited reserve the right to change or update this literature without prior notice.
Registered in: England and Wales. Registered No: 1289798



