

smartsigns®

Compact 300 Series

HUNTLEIGH

Flexible, reliable
and affordable



...performance for life

Simplicity, Reliability and Affordability

The SC300 is a new mains / battery operated compact spot check vital signs monitor; this powerful, flexible platform offers healthcare professionals a range of measurements which includes non-invasive blood pressure, pulse rate, pulse oximetry and a wireless tympanic temperature option.

Equipped with a long life 12 hour high capacity battery pack, the SC300 can be used in any healthcare setting.

Patient status is displayed on the 6" colour coded digital display enabling clinicians to make critical decisions at the point of delivery.

The SC300 is suitable for adult, paediatric and neonatal patients.





The SC300 sits comfortably in any environment, hospital ward, consulting room or clinic

Maximise Productivity

The attractive design with its easy to use user interface maximises the care givers' productivity, measurements are made at the touch of a button which means that users are able to spend more time caring for their patient rather than working through user documentation or set up options.

The wireless thermometer option delivers accurate measurements with all patient groups in seconds. The pre-heated probe tip enhances the patient experience and supports accurate measurements.

To support patient administration and effective record keeping, the large internal memory can store up to 50 sets of measurements.

The products' smooth surfaces and sealed keypad aids cleaning thereby reducing the risk of any cross infection. The integrated carry handle assists with transportation and portability.

The SC300's flexibility is enhanced through a choice of mounting options, users can select either a mobile stand with integrated accessory storage, a wall mount option also with a storage capability or an IV pole attachment where space may be limited.



Wall mount option with basket



Wall mount option with hooks

IV pole attachment

smartsigns®

Compact 300 Series

“As a caregiver, today’s healthcare provision is all about being effective and efficient, the SC300 with it’s flexibility meets these challenges...”



Technical Specification

Smartsigns Compact 300 Series	
Size	125 x 299 x 130mm
Weight	1.25Kg
Display type / size	LED, 100 x 120mm
Power voltage	100 – 240VAC
Power frequency	50 / 60Hz
Input current	0.15 – 0.3A
Battery type / capacity	Lithium Ion, 11.1V, 2200 mAh
Thermometer battery type / capacity	LR03 (AAA x 2), 1.5VDC
Patient groups	Adult, Paediatric & Neonate
NiBP	Oscillometric
SpO ₂	0% to 100%, 1% resolution
Temperature (option)	Tympanic, 34°C to 42.2°C (93.2°F to 107.6°F)

Smartsigns Compact 300 Series - Order Codes	
Order Code	Description
SC300	NiBP, Pulse and SpO ₂
SC300T	NiBP, Pulse, SpO ₂ and Temperature
SC300N	NiBP, Pulse and SpO ₂ (Nellcor)
SC300NT	NiBP, Pulse, SpO ₂ (Nellcor) and Temperature

Nellcor is a trademark of Covidien AG.

HUNTLEIGH *...performance for life*

Huntleigh Healthcare Ltd. - Diagnostic Products Division

35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom

T: +44 (0)29 20485885 F: +44 (0)29 20492520 E: sales@huntleigh-diagnostics.co.uk

W: www.huntleigh-diagnostics.com



Registered No: 942245 England & Wales. Registered Office:
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF
©Huntleigh Healthcare Limited 2015

AN ARJOHUNTLEIGH COMPANY, MEMBER OF GETINGE GROUP

® and ™ are trademarks of Huntleigh Technology Limited

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

AW:1001005-1

LIT 781310/EN-1