# seca 813 bt

# Digital flat scale with bluetooth interface and high capacity



- Solid steel frame, up to 440 lbs capacity.
- Non-slip rubber base for safe standing during measuring.
- Large LCD figures for easy reading.
- Step-off function for weighing without switching scales.
- Transmits weight via Bluetooth®.



### seca 813 bt:

# High capacity digital flat scale for individual patient use.

eca has designed the seca 813 bt with a bluetooth interface to transfer the Auto-HOLD value automatically to an external device. The scale offers a high capacity of up to 440 lbs and an extra wide platform. The tough rubber coating takes the heaviest challenges lightly and the flat construction makes stepping on very easy.



# A scale of sturdiness – with a touch of elegance.

Thanks to its solid steel frame, the seca 813 bt can withstand 440 lbs and features large LCD figures for easy reading. Furthermore, its extrawide surface and flat shape facilitate trouble-free weighing of heavyweight patients.

## Pleasant and comfortable weighing.

All platforms are not created equal. For this weighing platform seca has chosen a high-quality, two-component nubby material which is soft like rubber and pleasant to the touch. But don't be fooled – this material is also tough. It fulfils two requirements at the same time: pleasant weighing and an easy-to-clean surface that stands up to all commercially available cleaning and disinfecting products.

#### Paying attention to details.

Weight and girth are no problem for the seca 813 bt. You don't even have to bend over to switch it on, thanks to the unique step-off function. An issue with many scales is the display. Thanks to the 1.1 inch-high LCD figures, it is extremely easy to read the measurements. You see, a scale from seca makes life more simple by paying attention to even the slightest details.



High-quality, two-component nubby material feels as soft as rubber.



Big LCD figures. For easy reading.

## seca 813 bt

### Technical data

- Capacity: 440 lbs / 200 kg
- Division: 0.2 lbs / 100 g
- Dimensions (WxHxD):
  17 x 1.9 x 14.7" / 433 x 47 x 373 mm
- Weight: 6.3 lbs / 2.9 kg
- Functions: auto-HOLD, step-off, lbs/kg switch-over, tip-on automatic switch-on, automatic switch-off
- Parts to be delivered:
  4 AAA batteries,
  4 extension feet for carpets
- Transmits weight via Bluetooth®