

ho CryoProbe™ x+
Beauty Renewed



CryoProbe™ x+

EXTRA PLUS - The ultimate cryosurgery device for intensive use and for all types of benign lesions.



- 1-3mm applicator
- 2-6mm applicator
- 4-8mm applicator
- 7-20mm applicator
- 1-3mm long applicator
- 2-6mm long applicator



BETTER TREATMENTS

Effective, no-anesthetic removal of small or large benign lesions, in just a few seconds.

FASTER RECOVERY

Perfect results for complete removal of all benign skin lesions.

EXTREME PRECISION

Microfine jet for precise control.
Safe for use on delicate skin or near the eyes.

EXTREME POWER

N₂O at high pressure (50-55 bars) and ultralow temperature (-89°C/-128°F) for effective cell destruction.

BETTER BUSINESS

Efficient, practical and extremely cost effective. For a quick return on investment. Ideal for podiatrists and dermatologists.

Delivered with: 1x CryoProbe X+ pen • 1x 16g N₂O cartridge • 6 applicators. Order additional 8g or 16g cartridges separately.