

CryoPen[®]



Reinventing Cryosurgery





LEADING-EDGE CRYOSURGERY

Discover CryoProbe and Freezpen, a range of innovative, powerful, hand-held cryosurgery devices, designed to completely remove all skin lesions and tumors with a low cost per treatment and extremely quick procedures.

Treatment is very fast and can be easily integrated into daily practice.

The extremity of the device emits a microjet of nitrous oxide (N₂O) under high pressure (50-55 bars), achieving in-tissue temperature of -89°C for complete tissue necrosis in seconds.

No anesthesia, no pain

CryoProbe / Freezpen treatment causes no pain, no damage to healthy tissue and little or no scarring. And it requires minimal post-treatment care.

Easy to use

If you have no experience performing cryosurgery, all you need are your diagnostic skills and a steady hand. If you are experienced in using traditional cryosurgery tools, you'll be delighted by exceptional performance, millimeter precision and ease of use!



SEE



THE CRYOPEN IN ACTION !

Treat all skin lesions

Treat common pathologies
in just seconds:
Adenoma, Distichiasis, Epulis,
Eyelid adenoma, Meibomian gland
dysfunction, Sarcoid, Warts, Tumors
and more....

Easy to manage

Each CryoProbe and Freezpen
device is supplied with its own range
of applicators. Additional cartridges
can be ordered. Cartridges have a
2-year
shelf life and can be conveniently
stored at room temperature.



PERFECT REMOVAL OF SKIN LESIONS

CryoProbe and Freezpen delivers exceptional results for the complete removal of all skin lesions.



Sebaceous adenoma



Cyst



Sarcoid



Perianal adenoma



Eyelid tumor



Ear tumor



Meibomian Gland Dysfunction



Epulis



Tumor



Adenoma



CHOOSE THE CRYOPEN THAT'S RIGHT FOR YOU

The CryoPen series includes a range of pen-like cryosurgery devices, applicators and cartridges to meet every practitioner's needs.

Your choice of device depends on your usage.

We will be happy to advise you.



CryoPen® M

Suitable for less-frequent use or for smaller/superficial lesions.
Takes 8g N2O cartridges.

- 1-3mm applicator |
- 2-6mm applicator ||

CryoPen® X+

Suitable for intensive and varied use.
Available in 2 sizes, for 8g or 16g N2O cartridges.

- 1-3mm applicator |
- 2-6mm applicator ||
- 4-8mm applicator |||
- 7-20mm applicator ||||
- 1-3mm long applicator |
- 2-6mm long applicator ||

CryoPen® B+

Suitable for frequent use or for larger/deeper lesions.
Takes 23.5g N2O cartridges.

- 1-4mm applicator |
- 2-6mm applicator ||
- 4-8mm applicator |||



Cartridges

8g N₂O

1 box contain 24 cartridges,
1 carton contain 12 boxes.

16g N₂O

1 box contain 6 cartridges,
1 carton contain 12 boxes.



CryoPen® o+

Suitable as an all-purpose cryosurgery device for regular use.

Available in 2 sizes, for 8g or 16g N2O cartridges.

- 1-3mm applicator |
- 2-6mm applicator ||
- 4-8mm applicator |||
- 7-20mm applicator ||||

CryoPen® XP

Suitable for intensive and varied use.

Available in 2 sizes, for 8g or 16g N2O cartridges.

- 1-3mm applicator |
- 2-4mm applicator ||
- 3-5mm applicator |||
- 5-8mm applicator ||||

23.5g N₂O

1 box contain 6 cartridges,
1 carton contain 12 boxes.

5 YEARS WARRANTY

2 years for applicators!



MANUFACTURED BY



H&O EQUIPMENTS
Reinventing Cryosurgery

USA Office

115 Fairchild St, Suite 370,
Charleston, SC 29492
USA

Phone : +1 888-248-2838
cryoprobe@ho-equipments.com

Operating hours
Monday – Friday
9:00-17:00 (GMT -4)

Head Office

Rue des Journaliers 1
7822 Ghislenghien
Belgium

Phone : +32 68 26 86 00
info@ho-equipments.com

Operating hours
Monday – Friday
9:00-17:00 (GMT +2)

www.ho-equipments.com
info@ho-equipments.com

DISTRIBUTED BY



SIGMA MEDICAL

Dedicated to Health

Athens Offices

Head Offices

Ippokratous 142
19004 Spata
Phone: +30210 6630520
mail@medical.gr

Show Room

Ag. Thoma 22
11527 Athens
Phone: +302107473149
athens@medical.gr

Thessaloniki Offices

K. Melenikou 7
54635 Thessaloniki
Phone: +302310 500400
thessaloniki@medical.gr

<https://medical.gr/>