



AFD Technology Diaphragm

Adscope® 600GP
Gold-platedAdscope® 600ST
Tactical

Ordering Information

Item	Description
600**	Adscope® 600, Cardiology, (Specify Color)
600GP	Adscope® 600GP, Cardiology, Gold-plated (Black Tubing Only)
600ST	Adscope® 600, Cardiology, Tactical (Black Tubing Only)

Color Choices for 600** Only



NOTE: To ensure correct color when ordering, add color suffix to model number (e.g., 600**DG**)



Adscope® 600 Platinum Edition

AFD technology and a sculpted chestpiece make this the ultimate acoustic stethoscope.

The 600 series features

- ADC's adjustable frequency design (AFD) technology for the acoustic response of a traditional bell and diaphragm in a convenient one-sided chestpiece
- Light pressure for low, bell-type frequency, firm pressure for high, diaphragm-type frequency
- Sculpted chestpiece precisely CNC machined to exacting tolerances from surgical stainless steel for outstanding performance, rugged durability, and impressive ergonomics
- Non-chill diaphragm retaining rim for patient comfort
- Cardiology headset with large-bore stainless-steel binaurals, fixed at a 15° angle to maximize comfort and acoustic seal
- Bi-lumen (dual bore) 19" PVC tubing stays soft and flexible
- Adsoft™ Plus snap-on silicone eartips for the ultimate in wearing comfort and acoustic seal
- Tactical model uses black anodized aluminum aural tubes
- Accessory kit includes additional Adsoft™ Plus eartips and scope ID tag
- Weighs 7.2 oz (203 g); weight may vary depending on finish
- Overall length 27"
- Assembled, inspected, and packaged in the USA
- Individually acoustic tested in the USA
- Lifetime warranty and free parts forever