



AFD Technology Diaphragm

## Adscope® 601 Convertible

Combining AFD technology with a convertible design makes this the most versatile of our cardiology scopes.

### The 601 series features

- Convertible chestpiece for adult diaphragm/bell or adult/pediatric diaphragm configurations
- ADC's adjustable frequency design (AFD) technology provides enhanced acoustic response in both adult and pediatric diaphragm
- Stainless-steel chestpiece and aural tube construction (tactical model has aluminum binaurals)
- Extra-deep bell for unsurpassed low-frequency response
- Non-chill bell and diaphragm retaining rims for patient comfort
- Cardiology headset with large-bore stainless-steel binaurals fixed at a 15° angle to maximize comfort and acoustic seal
- Bi-lumen (dual bore) 19" PVC tubing stays soft and flexible
- Adsoft™ Plus snap-on silicone eartips for the ultimate in wearing comfort and acoustic seal
- Includes accessory kit with adult bell chestpiece, extra Adsoft™ Plus eartips, and scope ID tag
- Weighs 7.6 oz (215 g) with adult bell and 8.6 oz (243 g) with pediatric diaphragm; weight may vary depending on finish
- Overall length 28"
- Assembled, inspected, and packaged in the USA
- Individually acoustic tested in the USA
- Lifetime warranty and free parts forever



Convert pediatric diaphragm to adult bell in seconds.

### Ordering Information

Item	Description
<b>601**</b>	Adscope® 601, Cardiology Scope (Specify Color)

#### Color Choices



**NOTE:** To ensure correct color when ordering, add color suffix to model number (e.g., 601**DG**)

