



AFD Technology Diaphragm

Adscope® 608 Convertible

AFD technology and convertible design make this the most versatile of our clinician scopes.

The 608 series features

- Convertible chestpiece for adult diaphragm/bell or adult/pediatric diaphragm configurations
- Extra-large bell (1 3/4" diameter) for unsurpassed low-frequency response
- Ultra-sensitive diaphragm (1 1/4" diameter) for amplification and crisp high-frequency response
- AFD technology further enhances acoustic response: use light pressure for low frequencies, firm pressure for high frequencies
- Non-chill bell and diaphragm retaining rims provide patient comfort
- Clinician headset boasts reinforcing yoke molded into flexible 22" PVC and stainless-steel binaurals fixed at a 15° angle for greater comfort
- Adsoft™ Plus silicone eartips for the ultimate in wearing comfort and acoustic seal
- Accessory kit with spare Adsoft™ Plus eartips and scope ID tag
- Weighs 5.8 oz (164 g)
- Overall length 30"
- Assembled, inspected, and packaged in the USA
- Individually acoustic tested in the USA
- Lifetime warranty and free parts forever

Convert pediatric diaphragm to adult bell in seconds



Ordering Information

Item	Description
608**	Adscope® 608, Convertible Clinician Scope, (Specify Color)

Color Choices for 608** Only

BK	BD	COP	DG	FG	FL
FRB	FS	FV	G	IMCA	IMRS
IST	LV	MCB	MCA	MG	MRS
N	O	P	PBCA	R	RB
ST	TL	TQ	V		

NOTE: To ensure correct color when ordering, add color suffix to model number (e.g., 608**DG**)

