

Pediatric Adscope® 604

Legendary acoustics, performance, and functionality for pediatric care.

The 604 features

- Combination pediatric chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- Extra-large bell (1½" dia.) for unsurpassed low-frequency response
- Ultra-sensitive 1½" diameter diaphragm for greater amplification and crisper high-frequency transmission with adjustable frequency design (AFD) technology
- Non-chill bell and diaphragm retaining rims for patient comfort
- Reinforcing yoke molded into flexible 22" PVC "Y" tubing
- Stainless-steel binaurals fixed at 15° angle for greater comfort
- Adsoft™ Plus snap-on silicone eartips in small and large
- Spare eartips and scope ID tag included
- Weighs 4.6 oz, overall length 29.5"
- Inspected, assembled, and packaged in the USA
- Individually acoustic tested in the USA
- Lifetime warranty and free parts forever

Ordering Information

Item	Description
604*	Adscope® 604 Pediatric Scope (Specify Color)

Color Choices



NOTE: To ensure correct color when ordering, add color suffix to model number (e.g., 604**P**)



AFD Technology Diaphragm



Infant Adscope® 605

An Adscope® for the smallest patients.

The 605 features

- Combination infant chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- ¾" diameter bell provides unsurpassed low-frequency response
- Ultra-sensitive 1½" diameter diaphragm for greater amplification and crisper high-frequency transmission
- Non-chill bell and diaphragm retaining rims for patient comfort
- Reinforcing yoke molded into flexible 22" PVC "Y" tubing
- Adsoft™ Plus snap-on silicone eartips in small and large for the ultimate in wearing comfort and acoustic seal
- Extra diaphragm, spare eartips, and scope ID tag included
- Weighs 4.1 oz, overall length 29.5"
- Inspected, assembled, and packaged in the USA
- Individually acoustic tested in the USA
- Lifetime warranty and free parts forever

Ordering Information

Item	Description
605*	Adscope® 605 Infant Scope (Specify Color)

Color Choices



NOTE: To ensure correct color when ordering, add color suffix to model number (e.g., 605**G**)

