



## Adscope® 618 **Animals** Scope

Bring a smile to little faces while delivering the serious acoustics you need for assessments.



AFD Technology Diaphragm

### The 618 series features

- Distinctive, molded-resin animal face chestpiece optimally sized (1½") for pediatric use
- Adjustable frequency design (AFD) technology switches between low, bell-type frequency and high, diaphragm-type frequency simply by altering pressure
- Non-chill bell and diaphragm retaining rims for patient comfort
- Clinical headset boasts reinforcing yoke molded into the flexible 22" PVC tubing
- Lightweight aluminum binaurals in matching satin finish
- Adsoft™ Plus snap-on silicone eartips for the ultimate in wearing comfort and acoustic seal
- Seven interchangeable animal face snap-ons included:
  - Panda
  - Frog
  - Monkey
  - Koala
  - Deer
  - Tiger
  - Bear
- Accessory kit with spare eartips and a scope ID tag
- Weights: 3 oz (85 g)
- Overall length 30"
- Inspected, assembled, and packaged in the USA
- Individually acoustic tested in the USA
- Lifetime warranty and free parts forever

### Ordering Information

Item	Description
<b>618*</b>	Adscope® 618 Animals®, Pediatric (Specify Color)
<b>618KIT</b>	Adscope® 618 Animals®, Replacement snap-on faces

Color Choices



**NOTE:** To ensure correct color when ordering, add color suffix to model number (e.g., 618**SF**)

