



845-1 Series



845-2 Series

## Adcuff+ Cuffs

Our all-new line of one-piece, reusable blood pressure cuffs is designed for measurement accuracy, patient comfort, and practitioner ease.

- Bladderless design is simple to use and easy to clean
- Our most comprehensive Size Guide marking system features three different indicators (color coding, size icons, and cm range) for rapid identification of correct cuff size
- Convertible design for use with single- or double-tube systems
- Hook-and-loop adhesive rated for tens of thousands of cycles
- Gauge hang tab (on most sizes) compatible with conventional pocket aneroids
- Surface mounted tube ports permit a more uniform fit around patient limb
- Not made with natural rubber latex or phthalates
- Three-year warranty

## Adcuff™ Cuffs and Bladders 1 Tube / 2 Tube

- Durable nylon construction
- Exclusive Size Guide™ marking system
- 1-tube inflation bladder with 50cm tubing length (25cm on 2-tube)
- For use with most palm, traditional pocket, clock, or mercury aneroid instruments
- Adult size available in 17 color options
- Individually packaged
- Three-year warranty
- Luer connector (804CN) not included (1 Tube)
- Luer connector (891F) not included (2 Tube)

### Ordering Information

1 Tube Ordering Information					
Infant	Child	Sm. Adult	Adult	Lg. Adult	Thigh
845-7IBK-1	845-9CBK-1	845-10SABK-1	845-11ABK-1	845-12XBK-1	845-13TBK-1
845-7IN-1	845-9CN-1	845-10SAN-1	845-11AN-1	845-12XN-1	845-13TN-1
845-7IOR-1	845-9COR-1	845-10SAOR-1	845-11AOR-1	845-12XOR-1	845-13TOR-1
845-7IOR-1	845-9CGR-1	845-10SARB-1	845-11AN-1	845-12XBD-1	845-13TBR-1

2 Tube Ordering Information					
Infant	Child	Sm. Adult	Adult	Lg. Adult	Thigh
845-7IBK-2	845-9CBK-2	845-10SABK-2	845-11ABK-2	845-12XBK-2	845-13TBK-2
845-7IN-2	845-9CN-2	845-10SAN-2	845-11AN-2	845-12XN-2	845-13TN-2
845-7IOR-2	845-9COR-2	845-10SAOR-2	845-11AOR-2	845-12XOR-2	845-13TOR-2
845-7IOR-2	845-9CGR-2	845-10SARB-2	845-11AN-2	845-12XBD-2	845-13TBR-2

### Additional Adult Cuff Color Options



**NOTE:** To ensure correct cuff color when ordering, add color suffix to model number (e.g., 845-11ARB-2)



Available in the colors and sizes below.  
Also available as an option with select Diagnostix sphygs.

### Ordering Information

1 Tube Ordering Information					
Infant	Child	Sm. Adult	Adult	Lg. Adult	Thigh
8451-7IOR-1	8451-9CGR-1	8451-10SARB-1	8451-11AN-1	8451-12XBD-1	8451-13TBR-1

2 Tube Ordering Information					
Infant	Child	Sm. Adult	Adult	Lg. Adult	Thigh
8451-7IOR-2	8451-9CGR-2	8451-10SARB-2	8451-11AN-2	8451-12XBD-2	8451-13TBR-2