



e-sphyg 3™ Digital Blood Pressure Monitor

Our all-new, clinical-grade digital blood pressure monitor is the economical, easy-to-use solution you've been waiting for.

The e-sphyg 3 features

- Auto mode for quick readings in as little as 30 seconds
- Manual mode for auscultation with a stethoscope
- Customizable three-reading average mode with individual measurement display
- Last-reading memory
- Irregular heartbeat detection
- Oscillometric measurement with systolic and diastolic blood pressures, pulse rate, and mean arterial pressure (MAP)
- High-capacity rechargeable battery, standard
- For use with adult and pediatric patients
- Validated for pregnancy and preeclampsia
- Includes Small Adult, Adult, and Large Adult Adcuffs with exclusive Size Guide Marking System to ensure accurate readings on virtually any patient
- Desktop design is self-standing and features built-in handle; or choose optional desk, wall, or mobile platforms for added versatility in any practice
- Three-year warranty on unit, cuffs, and platforms; one-year warranty on battery and accessories

Ordering Information

Unit Models

e-sphyg 3 with Sm Adult, Adult, Lrg Adult Cuffs	9003K-MCC
e-sphyg 3 with Sm Adult, Adult, Lrg Adult Adcuff ⁺	9003K-MCC1

See page 45 for information on our new one-piece Adcuff⁺ line.

Mounting and Platforms

e-sphyg 3 Desk Caddy	9003D
e-sphyg 3 Desk Caddy for 2200 and 429	9003TS1
e-sphyg 3 Desk Caddy for 2200 and 433	9003TS2
e-sphyg 3 Mobile Stand and Basket	9003M
e-sphyg 3 Wall Mount and Basket	9003W
e-sphyg 3 Padded Carrying Case	9003BAG

Customers ordering in Canada can call 1-800-ADC-2670 for model availability.

Included Cuff Sizes:

Small Adult	Adult	Lg. Adult
-------------	-------	-----------



Features & Comparison

Feature	CLOCK ANEROID	E-Sphyg 3	Adviev 2
Accuracy ¹	+/- 3mmHg or less	+/-3mmHg or less	+/- 5mmHg or less
Accuracy Maintenance	Must periodically test against a reference standard	Self test lets operator know if calibration is needed ²	Self test lets operator know if calibration is needed ²
Display	Analog	LCD	LCD
Inflation	Manual	Automatic	Automatic
Deflation	Manual	Automatic	Automatic
Measurement Method	Auscultatory (Manual with scope)	Oscillometric (in Auto Mode) Auscultatory (in Manual Mode)	Oscillometric (in Auto Mode) Auscultatory (in Manual Mode)
Heart Rate	No	Yes	Yes
MAP	No	Yes	Yes ³
BP Averaging	N/A	Yes	Yes
Automatic Inflation	N/A	Yes (To Automatic or Preset Levels)	Yes
Manual Blood Pressure Mode	Yes	Yes	Yes
One-Button Inflation	No	Yes (To Preset Levels)	Yes (Fully Automatic)
Automatic Deflation	No	Yes	Yes (3-4 mmHg/sec)
Cuff Sizes Compatible with Auto Mode	N/A	(3) Sm. Adult, Adult, Lg. Adult	(5) Child, Sm. Adult, Adult, Lg. Adult, Thigh
Cuff Sizes Compatible with Manual Mode	6	3	5
Cuffs Included	Adult	3 Sm. Adult, Adult, Lg. Adult Traditional or Adcuff+	2 Adult, Lg. Adult Adcuff+
Power Supply	Manual	NiMh Battery/ AC Power	Li-ion/ AC Power
Readings on Single Charge	N/A	800-1200	200
Platforms			
Bag	No	Yes	No
Desk	Yes	Yes (2 styles)	Yes
Wall	Yes	Yes	Yes
Mobile	Yes	Yes	Yes
Memory	No	Last Reading	99 Reading
Upgradeability (Measures other parameters)	No	No	Yes Optional Temp & SpO2 Modules
Printer or PDF Interface	N/A	No	Yes
EMR Interface	No	No	EMR Capable
Calibration Warranty	Generally Lifetime	3 Year	3 Year

¹ Accuracy of pressure sensor

² Periodic check against a reference standard is recommended

³ Not validated for MAP in USA