

memmert

trust the best



Water bath

**WTB**

## WTB | PRACTICAL AND ROBUST

### Easy to clean

- rounded corners
- Stainless steel surfaces

### Easy to empty

- new drain system

### Easy to handle

- smallest possible footprint




## WTB | INTUITIVE AND USER-FRIENDLY

### Easy to operate

- 3.5 inch touchscreen
- New user interface

### Versatile functions



- Timer
- Alarm
- 2-point calibration
- Setting Celsius / Fahrenheit



memmert

## WTB | CUSTOMVIEW

Thanks to the customisable remote View **CustomView** you always keep an eye on the important information.



## WTB | TECHNICAL DATA



**WTB6**  
7.5 l



**WTB11**  
10 l



**WTB15**  
15 l



**WTB24**  
24 l



**WTB35**  
37.5 l



**WTB50**  
51 l

Model sizes / Article description	6	11	15	24	35	50
Power consumption at 230/115 V, 50/60 Hz (W)	1000		1420		2800	
Operating temperature range (°C)	min. 5 above ambient temperature up to +100					
Setting temperature range (°C)	+10 to +100					