## PRODUCT SPECIFICATION



# CEDERROTH

## 1 GENERAL PRODUCT DESCRIPTION

#### Description

A First Aid station containing a variety of first aid products, intended to be used in first aid situations, especially suitable for work places where there is a risk of burns. All products are easily available and protected against dust and moisture by transparent lids. Space for additional products is provided behind a lid with first aid instructions. Blue adhesive bandages (plasters) and blue soft foam bandage are placed in dispensers for easy access and theft-proof storage. A special key is provided for refill replacement. Medical devices in the station are CE-labelled.

The station contains >70% post-consumer recycled plastic (>85% in green coloured parts).

# **Storage information**

Dry storage in room temperature. For maximum protection, keep in original package until use.



#### 2 DETAILED PRODUCT INFORMATION

# Content per consumer unit

<u>Product</u>	<u>Qty</u>	<u>REF</u>	<u>Product</u>	<u>Qty</u>	<u>REF</u>
		(refill)			(refill)
First Aid station (empty)	1	-	Burn Gel Dressing 2-p	2	901900
Cederroth 4-in-1	2	1910	Burn Gel Spray	1	51011005
Cederroth 4-in-1 Mini	3	1911	Burn Cover	1	51901903
Salvequick Blue Detectable Plasters, regular, in dispenser	1	51030127	Cederroth Soft Foam Bandage, blue, 6cm x 2m	1	51011011
Salvequick Blue Detectable Plasters, regular/fingertip, in dispenser	1	51030126	Salvequick Wound Cleanser 20-p	1	323700
Cederroth Breathing mask	1	51616900	Key to the dispensers First Aid Instruction (label)	1 1	51616620 -

Dimensions (w x d x h): (station) 558 x 286 x 117 mm

Weight: (full station) 2943 g

Material

First aid station: Green parts: ABS-plastic, made of >85% post-consumer recycled ABS

Transparent lid: Polycarbonate

Sealing profiles: EPDM rubber (synthetic rubber)

Packaging: >90% of paper-based material, which is biobased and renewable

First aid products: Refer to Product specifications for each product

**Shelf life:** Refer to Product specifications for each product in the station.

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610026936	7310800015412	7312020016026
Content:	1 Station	1 CU	11 SKU/layer x 3 layers
Size (w x l x h mm):	294 x 572 x 123	294 x 572 x 123	800 x 1200 x 1032
Gross weight:	3268 g	3268 g	133 kg
Packing material:	Carton 313 g	-	-
	LDPE plastic 12 g		



Version: 03 Issued: 2022-02-04