

HEINE mini3000® F.O. Diagnostic Sets



Set complete with: **mini3000 F.O. Otoscope**
 1 set (4 pcs.) reusable tips (B-000.11.111)
 5 each of 2.5 and 4 mm dia. AllSpec disposable tips
 hard case

mini3000 F.O. Diagnostic Sets	2.5V XHL
with mini3000 battery handle with batteries	D-851.10.021



Set complete with: **mini3000 Ophthalmoscope, mini3000 F.O. Otoscope,**
 1 set (4 pcs.) reusable tips (B-000.11.111)
 5 each of 2.5 and 4 mm dia. AllSpec disposable tips
 hard case

mini3000 F.O. Diagnostic Sets	2.5V XHL
with 2 mini3000 battery handles with batteries	D-873.11.021




Set complete with: **mini3000 Ophthalmoscope, mini3000 F.O. Otoscope,**
 mini NT table charger
 5 each of 2.5 and 4 mm dia. AllSpec disposable tips

mini3000 F.O. Diagnostic Sets	2.5V XHL
with 2 mini3000 rechargeable handles	D-859.11.022*

*with this Art. No. only available in black.

HEINE AllSpec® Disposable Tips



- ⊗ **For the most secure fit, best light transmission and hygiene, we recommend the use of HEINE Otoscope Tips.**
- ⊗ **Assured quality.** No sharp edges.



AllSpec Disposable tips (grey) for BETA400, BETA200, K180, mini3000 F.O., mini3000 Oscopes	
Pack of 1000 pcs. 2.5 mm dia. (infants) [03]	B-000.11.128
Pack of 1000 pcs. 4 mm dia. (adults) [03]	B-000.11.127
Case of 10 dispenser packs of 250 tips, 2.5 mm dia (infants) [01]	B-000.11.151
Case of 10 dispenser packs of 250 tips, 4 mm dia (adults) [01]	B-000.11.150
AllSpec Tip dispenser, white [02]	B-000.11.149

- ⊗ **AllSpec Tip dispenser** incl. tips Ø 2.5mm and 4mm [02].
The right size, visible and instantly available. A clean new tip for each and every patient.

⊗: Single use only.



Reusable Tips




Ø 2,4 mm Ø 3 mm Ø 4 mm Ø 5 mm

- ∴ For the most secure fit, best light transmission and hygiene, we recommend the use of **HEINE Otoscope Tips**.
- ∴ **Made of high-density SANALONS plastic with smooth external finish.** Easy to clean. All common germicides can be used. Autoclavable at up to 134 °C.
- ∴ **Ergonomic shape.** Painless atraumatic insertion.

Reusable tips (black)	1 piece	50 pieces
for BETA 400, BETA 200, K180, mini3000 F.O., mini3000 Otoscopes		
Reusable tip 2.4 mm dia.	B-000.11.107	B-000.11.207
Reusable tip 3.0 mm dia.	B-000.11.108	B-000.11.208
Reusable tip 4.0 mm dia.	B-000.11.109	B-000.11.209
Reusable tip 5.0 mm dia.	B-000.11.110	B-000.11.210

1 set = 4 pcs. reusable tips SANALONS 2.4/3/4/5mm dia. **B-000.11.111**

Soft Disposable Tips



Ø 3 mm Ø 5 mm

- ∴ Fits over the end of the tip to protect the auditory canal.

Soft disposable tips (black)	
Soft disposable tip 3 mm dia. 40 pcs./1 pack	B-000.11.141
Soft disposable tip 5 mm dia. 40 pcs./1 pack	B-000.11.142

⊗: Single use only.

Insufflation Bulb



- ∴ For pneumatic testing.

Insufflation bulb for BETA 400, BETA 200 and K180 Otoscopes	B-000.11.240
Insufflation bulb for mini3000 F.O. Otoscope with connector	D-000.80.102
Insufflation port connector for mini3000 F.O. Otoscope	D-000.80.101