

# MIO-SONIC

Mio-Sonic is the compact and powerful ultrasound device, very suitable for muscle contractures and joint stiffness.

## MEDICAL DEVICE CLASS IIa

MIO-SONIC is a medical device CE0068

## KIT

- ✓ Mio-Sonic
- ✓ Medical power supply
- ✓ Ultrasound gel
- ✓ User manual
- ✓ Carriage and storage bag
- ✓ Box

## TECHNICAL FEATURES

- ✓ Professional device for ultrasoundtherapy
- ✓ 9 programs
- ✓ 1 MHz ultrasound, maximum output power 6.4W, maximum intensity 1,6W/cm<sup>2</sup>
- ✓ Collimated ultrasound beam for an optimal penetration beam
- ✓ 3 ultrasound beam intensity levels: 0,08 W/cm<sup>2</sup> (low) – 0,8 W/cm<sup>2</sup> (medium) – 1,6 W/cm<sup>2</sup> (high)
- ✓ Medical power supply
- ✓ Class IIa medical device
- ✓ Direct contact device (not suitable for water immersion).

## DIMENSIONS

200 x 47 x 69 mm

## WEIGHT

190 g



# MIO-SONIC

## PROGRAMS

### TREATMENTS

Osteoporosis	Tendinitis	Bursitis	Drainage
Arthritis	Pubalgy	Brachialgy	Modelling
Neuralgias	Tallonitis	Capsulitis	Stretch marks
Cervicalgias	Contusions	Chondropathy	Microlifting
Sciatalgy	Sprains	Cruralgy	Wrinkles
Lumbalgy	Hematoma	Acne	Phonophoresis
Carpal Tunnel	Oedema	Venous stagnation	Arms cavitation
Epycondilitis	Lymphedema	Muscular contractures	Abdomen cavitation
Epitrochleitis	Adhesions	Cramps	Limb cavitation
Periarthritis	Algodystrophy	Massage	Gluteus cavitation

## HIGH PERFORMANCE DEVICE

Innovative and high performance device both for home use and professional use. It combines compactness and high performances in a small, light and easy to use applicator.

The lightness and handiness allows a continuous use without any difficulty.

Mio Sonic can be used for the following applications:

- ✓ Pain therapy
- ✓ Fitness treatments
- ✓ Esthetic treatments