



HEINE DELTA® 20 T Dermatoscope

Highest Quality Images and Diagnostic Flexibility – with LED^{HQ}



LED NOW IN HEINE QUALITY.



EXCLUSIVELY AT HEINE

- Toggle function
- High performance LEDs
- High-resolution optics

"Blink Sign" observable in a basal cell carcinoma



Basal cell carcinoma non-polarised illumination



Basal cell carcinoma polarised illumination

For dermatoscopic examination for early diagnosis of malignant melanoma. Also for the examination of non-melanocytic lesions, basal cell carcinoma and dermatofibroma.

- **Precision HEINE Optics.** High quality, achromatic optical system provides unmatched sharpness and resolution. Focusing system provides for 10x to 16x magnification for recognition of even the smallest details and structures.
- **Adjustable Eyepiece.** Individual focusing, range of correction approx. -6 to +6 dpt.
- **LED in HQ:** The new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for the most accurate diagnosis. Red is red, blue is blue. Colour temperature: typ. 5,000K, Colour Rendering Index (CRI) >80.
- **Optional small contact plate available.** Special small contact plate for difficult to access lesions.
- **Digital documentation capabilities.** Connect to major brand SLR digital cameras with HEINE digital camera lens system. Specially designed lenses integrate the camera and DELTA20 T optics together to provide the highest quality imaging possible.
- **Toggle between polarised and non-polarised illumination modes.** Conveniently alternate between polarised and non-polarised views with the press of a button – no exchange of contact plates necessary. Allows recognition of the "Blink Sign" when observing crystalline structures and presence of milia cysts, increasing diagnostic function.

DELTA20 T Dermatoscope Head	LED
with contact plate Ø 23mm with scale and dermatoscopy compendium, without handle	K-008.34.221
with contact plate Ø 23mm without scale and dermatoscopy compendium, without handle	K-008.34.222

HEINE DELTA® 20 T Dermatoscope Sets



LED NOW IN HEINE QUALITY.



Set complete with: **DELTA20 T Dermatoscope**
contact plate 23mm Ø dia. with scale
10ml dermatoscopy oil
dermatoscopy compendium
hard case

DELTA20 T Sets	LED
BETA battery handle (2.5V)	K-262.10.118
BETA4 USB rechargeable handle with USB cord and plug-in power supply	K-262.28.388
BETA4 USB rechargeable handle	K-262.28.387
BETA4 NT rechargeable handle with NT4 table charger	K-262.24.420

Accessories

For DELTA 20 T Dermatoscope



Contact plates for DELTA 20 T

23mm dia. with scale [01]	K-000.34.216
23mm dia. without scale [01]	K-000.34.217
Small contact plate, 8 mm dia. [02] for the examination of inaccessible lesions	K-000.34.205



Dermatology oil

(pack of 6 x 10ml). Reduces reflexes on the outer layers of skin.	K-000.34.005
---	--------------



SLR Photo adaptor

for Dermatoscope/digital SLR camera, with optics.
Wide Field of View, crisp image

for Canon	K-000.34.185
for Nikon	K-000.34.186



Adaptor cord

(1.3 m) for instrument/handle	X-000.99.231
-------------------------------	--------------



BETA Belt Clip

for BETA handle	X-000.99.008
-----------------	--------------



Photo Accessory Set

includes SLR Photo adaptor, adaptor cord, distance ring and BETA Belt Clip

for Canon	K-000.34.190
for Nikon	K-000.34.191