



Dr. pen[®] A9






MEDICAL.GR
Dedicated to Health

LED LIGHT GEAR SPEED DISPLAY 6 SPEEDS CONTROL

Speed control system, increase speed control, increase speed.



MEDICAL.GR
Dedicated to Health


SPECIAL PACKAGING

Convenient storage. Each accessory can be hold with exclusive place.



MEDICAL.GR
Dedicated to Health

NEWLY UPGRADED ELECTRIC MICRONEEDLE



LED light gear speed display
6 speeds control

Speed control system,
increase speed control,
increase speed.



Ultra-long working time

Big capacity battery to ensure the
pen working efficiently.



Dr. pen® A9



Special packaging

Convenient storage. Each accessory can be hold
with exclusive place.

Valve needle cartridge

New design needle cartridge. Working efficient and
stable with 3 bayonet lock design. Build-in safety
isolation zone to ensure the sealing performance.
More hygiene and safe.



MEDICAL.GR
Dedicated to Health



VALVE NEEDLE CARTRIDGE

New design needle cartridge. Working efficient and stable with 3 bayonet lock design. Build-in safety isolation zone to ensure the sealing performance. More hygiene and safe.



MEDICAL.GR
Dedicated to Health



MEDICAL.GR
Dedicated to Health



ULTRA-LONG WORKING TIME

Big capacity battery to ensure the pen working efficiently.



MEDICAL.GR
Dedicated to Health





Model	Dr.pen A9	Packing List
Power Model	Wirless/Wired Model	① Box ② Tray ③ Dr.Pen A9 ④ protective cover ⑤ Manual ⑥ Cartridges *2 pcs ⑦ Adapter ⑧ Cable
Speed Level	6 Speeds	
Needle Adjustment Range	0-2.5mm	
Adapter	US/EU/UK/AU/KR	
Handle Material	ABS	
Package Size	169*123*56mm	
Package Weight	420g	



MEDICAL.GR
Dedicated to Health

Dr. pen® A9



Specification

Brand	Dr.pen
Model	A9
Power model	Wireless & Wired
Battery capacity	700 mAh (wireless model)
Speed	6 levels
Adjustable needle length	0 - 2.5 mm
Pen material	ABS
Available adapter	US/ EU/ UK/ AU/ KR
Input	AC 100-240V, 50/ 60Hz
Output	DC 5V, 1A
G.W.	17*13*6 cm/ 1 pc/ 420 g
1 carton	38*29*35 cm/ 20 pcs/ 9 kg
Customization (OEM & ODM)	logo, private label, manual, package, new shape, change color.

Advantage



- LED light gear speed display with incremental speed control mode.
- New design Valve needle cartridge, built-in safety isolation zone to ensure the sealing performance, more hygienic and safe.
- New appearance, minimalist and modern. Double colored plastic handle, wear resistant with comfortable grab feeling.
- Wireless model. Ultra-long working time, big capacity battery to ensure the pen working efficiently.
- New cartridge 18 needles, diameter 0.18mm. Less pain and better efficacy.
- Nice white travel box, convenient storage and portability. Each accessory can be held in exclusive space. Minimalist and practical.
- 1-6 speed levels can meet different needs of treatment.
- 3 bayonets lock design for needle cartridge, works efficiently and stably.

Packing list



Pen: 1 pc Cable: 1 pc Manual: 1 pc
 Color box: 1 pc Adapter: 1 pc White inner box: 1 pc
 Needle cartridge: 12 pins * 2 pcs

Needle cartridge for MTS

