











Needle Cartridge	Needle Cartridge Type	Needle Cartridge Pins	Features	Matched Device
A1	PMU 	R1/F3/R3/F5/R5/F7/F9	Blue, 2-point snap positioning	 Dr.pen A1
	MTS 	9/12/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R		
N2 Or M5 Or M7 Or E30	PMU 	R1/F3/R3/F5/R5/F7/F9	Transparent, 2-point snap positioning	 Dr.pen N2, M5, M7, E30
	MTS 	9/12/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R		
A6	PMU 	R1/F3/R3/F5/R5/F7/F9	Blue, 2-point snap positioning	 Dr.pen A6
	MTS 	9/12/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R		
A7	PMU 	————	Gray, Hidden Buckle	 Dr.pen A7
	MTS 	9/12/24/36/42/Nano-S/ Nano-3D/Nano-R		
X5	PMU 	————	light Blue, Screw Bayonet	 Dr.pen X5
	MTS 	9/12/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R		
M8	PMU 	————	Gray Cyan, 3-point snap positioning	 Dr.pen M8
	MTS 	11/16/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R		
A6s	PMU 	————	Green, 3-point snap positioning	 Dr.pen A6s
	MTS 	11/16/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R		
H2	PMU 	————	Transparent, Press to install the needle	 Hydra pen H2
	MTS 	12 pins:0.25mm/0.5mm/1mm Nano-HS/Nano-HR		
V6 Or V8 Or V9 Or V11	PMU 	M1(0.3mm)/L1(0.4mm) /F3/R3/F5/R5/F7/F9	White, Screw Bayonet	 Artemex V6 V8 V9 V11
	MTS 	9/12/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R	Blue	
V7	PMU 	M1(0.3mm)/L1(0.4mm) /F3/R3/F5/R5/F7	Transparent, Screw Bayonet	 Artemex V7
	MTS 	9/12/24/36/42/Nano-S/ Nano-3D/Nano-5D/Nano-R		