



### Automatic application of serum with adjustable liquid output

Hydra Pen uses essence liquid to penetrate the skin automatically when working.




Hydra Pen is an innovation technology in the field of aesthetic, it's also a new generation of MTS brand founded by Ekai in 2018, HydraPen takes not only the advantages from dermaroller and dermapen, more importantly, which works with its patent cartridge to apply the serum or essence automatically.

## Advantages

### Adjustable liquid output

In order to adapt to essence with different consistency, there are three levels of customized liquid output, and the liquid output can be controlled with the speed gear. The higher the speed gear, the more liquid output frequency, and vice versa.



### Three kind of needle types solve different skin problems

\* Needles are disposable, safe and sanitary



H12

Removing acne, Acne pits and marks, scar and stretch marks




Nano-HR

Whiten the skin, moisturize and tighten, shrink pores, and face fine lines



Nano-HS

For fine lines around eyes, sensitive muscles, red blood




\* Continuous use can improve the skin, and the actual effect varies from person to person



### Automatic application of serum

When the cartridge acts on the skin up and down, the microneedle pierces the epidermis and opens thousands of fine channels. At this time, the essence liquid will automatically penetrate from the applicator and reach the deep skin evenly, with the absorption rate up to 90% higher than that of manual application. With the function of microneedle, the skin tissue repair can be quickly started, stimulate cell regeneration, restore the vitality and elasticity of the skin, and double the healing effect.



### LED Screen Visual Digital Display



Due to the effect of force, the higher the speed gear, the more the frequency of liquid discharge, and vice versa.

### Multi-faceted care

Meet different skin care needs

Adjustable needle length 0-1.0 mm  
Turn the needle counterclockwise from short to long




H3 with low battery warning helps users to know the battery remaining power intuitively, Better control the working hours.

The number of speed flashes slowly when the battery capacity has 10% left

The pen will turn off automatically after the number of speed flashes quickly for 1 minutes.

The pen shows "L" when restarting after the low battery automatic shutdown.



## Product Parameters



### Specification

Brand	Hydra pen	Model	H3
Loadable liquid capacity	3 ml	Speed Level	6 Speeds
Adjustable liquid output	3 levels	Battery capacity	700mAh
Needle Adjustment Range	0-1.0 mm	Adapter	US/EU/UK/AU/KR
Handle Material	Metal	G.W.	17*12.3*5.6 cm/1 pc/ 390 g
Available Cartridges	H12, Nano HS, Nano HR	1 carton	38*29*35 cm/20 pcs/9 kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box      ② Hydra pen H3 ③ 12 pins * 2 pcs    ④ Adapter    ⑤ Cable ⑥ Manual    ⑦ Retaining ring* 2 pcs

