



DRS<sup>®</sup>  
Dermaroller  
System

Founded in 2010, DRS is the first microneedling brand based on the MTS mesotherapy, which break the traditional skin care theory and increase absorption 90%. DRS provides various models for different skin parts and different skin symptoms. In the past decade, they are loved by the customers from all over the world. Let us happy rolling!



## Advantages



- High quality 192 titanium alloy needle material manufacturing ensuring maximum durability, not easy to bend or break.
- Sharp and fine needles, penetrate the skin as less force, avoid skin discomfort.
- Seamless roller head, no hiding dirt, easy to disinfect and clean, more hygiene and safety.
- Protective storage case keeps needles safe and clean.
- Beneficial for hair and beard growth treatment.
- Unique and nice round bottle packaging, protected products well.

## Packing details

GTSROLLER 192 : 1 pc

Round bottled: 1 pc

Color box: 1 pc



## Specification

Brand	GTSROLLER	Model	192
Needle Number	192 pins	Needle length for choice	0.2 / 0.25 / 0.3 / 0.5 / 0.75 / 1.0 / 1.5 / 2.0 / 2.5 / 3.0 mm
Needle material	Titanium	G.W.	16*7*7 cm/ 1 pc/ 100 g
Handle/roller material	PC	1 carton	68*36*38 cm/ 100 pcs/ 11 kg
Color	Gold and silver	Customization (OEM & ODM)	logo, private label, manual, package, new shape