



bio roller® G5

Initiate skin repair Rejuvenate collagen skin

EMS micro current
540pcs micro needles for skin cell repair
Two LED lights with vibration care



Advantages



**EMS electroporation
micro needle**
Stimulate cells and accelerate
collagen synthesis

uniform conductivity

The microneedles penetrate the epidermis and open up thousands of tiny channels. The 9 needles continuously release weak currents to generate gentle nanoprop to penetrate deeply into the muscle layer, promote deep muscle movement, and reshape firm contours. Initiate skin tissue repair, stimulate cell regeneration, repair skin, shrink pores, and restore skin elasticity.

Two LED lights + Vibration care

After microneedling treatment, massage and soothe the skin with two different lights and vibrations to improve microcirculation.


- Blue light (455 nm)**
Calm and Soothe
- Red light (630 nm)**
Fade wrinkles and rejuvenate skin



Adjustable 5 Levels of EMS microcurrent



① Blister Packaging



③ Easy to install and remove



② Clean & sanitary



④ Face the skin Sign



Skin Care

Lift & Tighten Skin

Microneedling: Accelerate pigment metabolism, whiten and rejuvenate skin, remove yellow, remove acne marks, lighten melasma, reduce acne pits, improve stretch marks, shrink pores

EMS microcurrent: lighten fine lines, firm skin, boost collagen, lift skin, face shaping

Blue light + vibration: anti-inflammatory blemish, enhance beauty effect, shrink pores, soothe and calm

Red light + vibration: concentrated anti-aging, reduce wrinkles, accelerate microcirculation, rejuvenate skin and lighten lines



Lift and firm the skin




Fade acne scars



Lighten fine lines

Product Parameters



Specification

Brand	Bio Roller	Model	G5
Power Model	Wireless	Battery capacity	500 mAh
EMS microcurrent	5 levels	Adapter	US/EU/UK/AU/KR
LED light	Blue light, Red light	Roller needle length	0.25mm, 0.5mm, 1.0mm, 1.5mm, 2.0mm
Vibration	2 levels	Charging mode	USB Type-C interface
Handle Material	ABS	G.W.	17.3*10.3*4.2 cm/1 pc/ 350 g
Roller head	0.5 mm* 2 Pcs	1 carton	36 *27 *33 cm/20 pcs/8.5 kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color	Packing List	① Box ② Bio Roller G5 ③ Manual ④ 0.5 mm roller head * 2 Pcs ⑤ Adapter ⑥ Cable

