



TECHNICAL DATA

iBEAUTY
Mate

Technical data

General

Power consumption	150 W
Voltage	220-230 V
Frequency	50-60 hz
Display	7" LCD, touch-sensitive IPS 1024 x 600 px
Ventilation method	Air fan
Base dimensions	38,5 x 21 x 21,5 (cm)
Weight	7,4 kg
Standards / Certification	CE

Nano-needle mesotherapy

Cartridge type	with 0,1 mm nano-needles
Frequency	level 1: 25-30 insertions/sec
	level 2: 35-40 insertions/sec
	level 3: 45-50 insertions/sec
	level 4: 55-60 insertions/sec
	level 5: 65-70 insertions/sec
Line feed depth	0 - 3 mm
Head dimensions	15 x 2 cm
Head weight	77 g
Head material	aluminium
Cable length	1,8 m

Diamond microdermabrasion

Pump capacity	12-15 l/min
Pump pressure	0-2,2 bar
Types of abrasive tips	7 mm (D-240)
	8 mm (D-200)
	14 mm (D-180)
	14 mm (D-180)
	14 mm (D-120)
	17 mm (D-240)
	17 mm (D-180)
	20 mm (D-180)
	20 mm (D-75) FOR FEET ONLY
	22 mm (D-120)
Head dimensions	1,5 x 11 cm
Filtration type	pre-filter, liquid filter on the hose
Hose length	1,8 m

Oxybrasion

Pump capacity	25-28 l/min
Pump pressure	0-12 bar
Output maximum pressure	2 bar
Head dimensions	15 x 7 x 3 cm
Head weight	110 g
Tank capacity	7 cm ³
Head material	aluminium
Hose length	1,8 m

Oxygen infusion

Pump capacity	25-28 l/min
Pump pressure	0-12 bar
Output maximum pressure	2 bar
Head dimensions	15 x 7 x 3 cm
Head weight	110 g
Tank capacity	7 cm ³
Head material	aluminium
Hose length	1,8 m

LED

Number of diodes	7
Current for a single diode	18 mA
Colours	white, blue, green, red, orange, yellow, purple
Dimensions	3,7 x 19 cm
Head weight	130 g
Head material	Poliaceta
Cable length	1,8 m

Shallow kinetic massage


Pump capacity	25-28 l/min
Output maximum pressure	2 bar
Pump pressure	0-12 bar
Head dimensions	3,7 x 19 cm
Head weight	650 g
Head material	Poliacetal, stainless steel, aluminium
Filtration type	condenser on the hose
Hose length	1,8 m

Endermomassage

Pump capacity	12-15 l/min
Pump pressure	0-2,2 bar
Under-pressure	0- 53 cmHg
Large head diameter	12 cm
Large head weight	700 g
Small head diameter	7,5 cm
Small head weight	280 g
Head material	Poliacetal
Roller material	Stainless steel
Filtration type	pre-filter + liquid filter on the hose
Hose length	1,8 m




Base dimensions




21,5 cm

21 cm



38,5 cm

Contents of the set



- | | |
|--|--|
| 1. iBeautyMate control device | 8. Diamond tips for Microdermabrasion (10 pcs.) with filters |
| 2. Head for Oxybrasion and Oxygen infusion (aerograph) | 9. Line for Endermomassage with a filter (grey) |
| 3. Flexible pressure line for Oxybrasion and Oxygen infusion | 10. Two heads for Endermomassage with filters (5 pcs. each) |
| 4. Head for Shallow kinetic massage | 11. Pen for Nano-needle mesotherapy with a line, cartridge (2 pcs. in a set) |
| 5. Line for Shallow kinetic massage with a condenser | 12. LED head with an integrated cord |
| 6. Line for Microdermabrasion with a filter (transparent) | 13. Power supply cord |
| 7. Head for Microdermabrasion (3 pcs.) | 14. Organiser for heads |

iBEAUTY Mate

IBEAUTY MATE SPÓŁKA Z OGRANICZONĄ ODPOWIEDZIALNOŚCIĄ
(SHARE CAPITAL 20,000pln Polish Zloty)
Wrocławska 35-37
62-800 Kalisz
VAT: PL5272944522
REGON: 387726659
KRS: 0000874265