



AVAILABLE
NOW!

**The new
HEINE visionPRO HyMac 3
hyperangulated blade.**

New and different but always familiar.

Every now and again, you encounter airway situations that require you to look 'around the corner' where it is beneficial to use a slightly more curved blade. The new HEINE visionPRO HyMac 3 single-use blade is now available specifically for these special situations.

It is curved in the distal area, causing it to tilt sharply behind the base of the tongue, thus helping to lift the epiglottis 'almost by itself'.

- ▶ The visionPRO HyMac 3 can be used like an ordinary visionPRO Mac blade – therefore it can be used intuitively without any need for extra training
- ▶ A standard stylet is all that is needed to navigate the ET tube safely to its destination
- ▶ The complexity of the intubation process is not increased unnecessarily
- ▶ Fast & intuitive

THE FIRST HYPERANGULATED BLADE MADE OF UPCYCLED MATERIAL.

- ▶ Environmental protection due to CO₂ savings
- ▶ Sustainable production on an energy-saving machine with green electricity
- ▶ Reduction in material consumption:
Minimalist design & injection moulding process that produces very little waste

SAFE AND HYGIENIC.

- ▶ Clearly visible labelling
- ▶ Secure fitting, even in time-critical situations
- ▶ The blade design protects the camera arm from contamination, and enables fast and effective reprocessing
- ▶ Each blade is packed in a hygienically clean manner in its own individual packaging

visionPRO HyMac 3 single-use blades, box of 10

F-000.22.246

HEINE Optotechnik GmbH & Co. KG

Dornierstr. 6, 82205 Gilching, Germany
www.heine.com



We reserve the right to change specification without notice.

07/23. A-000.00.239 EN

