

## Hydrosun® 750 / 750N

The new standard in the care of prematures, newborns and adults!

Radiant heaters are used everywhere in the clinical care of premature and newborn babies. The focus is on the continuous maintenance of the core body temperature. Immediately available warmth - intense and gentle for the child - nature shows us! wIRA radiation is our concept: The surface-dangerous radiation areas that cause burns and allow the skin to dry out get filtered out, thanks to a patented water filter from the radiation spectrum.

Areas of application - guarantee of maximum flexibility. The Hydrosun® 750/750 N can be used to treat newborns on the changing table, as a therapy radiator used in the area of premature babies or in anesthesia. Despite the small radiation field, more heat can be applied to this heater without the risk of burns than with conventional ceramic heaters. Due to the small radiation field, the system can be used universally - whether to maintain the Body core temperature during examinations and interventions or as a red light emitter for skin diseases. The Treatments and the environment as well as devices used on children (e.g. incubators) are not, however heated.


- Above the changing table / examination area with wall or ceiling arm
- In functional diagnostics / X-rays
- Mobile for interventions / examinations
- In the delivery room for mother and child
- In anesthesia

### Benefits:


- **Full performance immediately**  
No preheating, no waste of energy, ready for immediate use
- **Simple, safe and effective technology**  
No alarms necessary
- **Gentle deep heating of the child**  
No warming up of the practitioners




Hydrosun® 750 N radiant heater



Heavy surface stress



Lower surface load



Great depth effect and rapid warming

### Ordering Information

Product Description	Ord.-No.
Radiant Heater Hydrosun® 750T (tabletop unit) without timer and filter, with spacer bar and 375watts spare bulb	HS750101T
Radiant Heater Hydrosun® 750N WH (for wall mount) without timer and filter, with spacer bar and 375 watts spare bulb	HS750101WH
Radiant Heater Hydrosun® 750WS (for wall bracket rotatable) without timer and filter, with spacer bar and 375 watts spare bulb	HS750101WS
Radiant Heater Hydrosun® 750NFS (for FS tripod) without timer and filter, with spacer bar and 375 watts spare bulb	HS750101FS
Radiant Heater Hydrosun® 750N PP (for PP tripod) without timer and filter, with spacer bar and 375 watts spare bulb	HS750101PP
FS tripod for Hydrosun® lamps of the 750 series (for wheelchairs/chairs)	HS750310
PP tripod for Hydrosun® lamps of the 750 series (for non-accessible beds/chairs)	HS750311
Wall arm, swivel for Hydrosun® 750 series systems*	HS750139
Ceiling arm, swivel for Hydrosun® 750 series systems*	HS750138
Black filter for Hydrosun® 300/305/500/750	HS0007
Damping filter for Hydrosun® Series 750/750N	HS750300
Hydrocuvette (replacement price) for Hydrosun® Series 750	HS750350
Halogen lamp for Hydrosun® 305N/750N incl. 2 fuses	HS0305

\* Delivery time about 12 weeks