

# SECTION DIRECT OPHTHALMOSCOPES

#### HEINE K180<sup>®</sup> Ophthalmoscope The standard direct ophthalmoscope A full featured instrument at an economical price. : Precision HEINE optics. Brilliant retinal image. : XHL Xenon Halogen Technology. Bright, white light. Clear retinal image. : Ergonomical shape. The instrument fits the orbita comfortably. : Soft orbital rest. Steadies the instrument and protects your glasses. : Choice of 5 apertures. For both large and small pupils. Choice of fixation star or cobalt blue filter for corneal examination. : Dust-resistant housing. Maintenance-free. : 28 lenses from - 35 D to + 40 D. Brilliant image even at high corrections. 1–10 15 20 40 + in 1 D steps 1-10 15 20 25 35 D - in 1 D steps : High-impact polycarbonate housing for durability. K180 Ophthalmoscope 2.5 V 3.5 V K180 Ophthalmoscope with aperture wheel 1 C-001.30.206 C-002.30.206 K180 Ophthalmoscope with aperture wheel 2 with cobalt blue filter C-001.30.207 C-002.30.207 XHL Xenon Halogen spare bulb X-001.88.084 X-002.88.086 · Reduces corneal reflex Dust-resistant

#### K180 apertures





Aperture wheel 2 with blue filter

## HEINE K180<sup>®</sup> Ophthalmological Diagnostic Sets







Set complete with K180 Ophthalmoscope, in standard version with aperture wheel 1 in soft pouch with:

	2.5 V	3.5 V
BETA battery handle	C-182.10.118	
BETA R rechargeable handle		C-182.20.376
BETA TR rechargeable handle and transformer		C-182.20.384
BETA NT rechargeable handle and NT 200 charger		C-182.20.420

### Set complete: K180 Ophthalmoscope, K180 Otoscope and 1 set of reusable tips in zipper case with:

	2.5 V
BETA battery handle	A-177.10.118

Set complete: K180 Ophthalmoscope, K180 Otoscope and 1 set of reusable tips in hard case with:

	2.5 V	3.5 V
BETA battery handle	A-279.10.118	
BETA R rechargeable handle		A-279.20.376
BETA TR rechargeable handle and transformer		A-279.20.384
BETA NT rechargeable handle and NT 200 charger		A-279.20.420