

ORIGINAL CU, PRIMEDIC AND PROGETTI PADS


Pair	5 pairs	PEDIATRIC PADS: BRAND - MODEL	Pair	10 pairs	ADULT PADS: BRAND - MODEL
33569	-	CU Medical System CU-ER1, CU-ER2, CU-ER3	35345	33565	CU Medical System CU-ER1, CU-ER2, CU-ER3
33443	-	iPad NF 1200	33441	33441	iPad NF 1200
35345	-	CU-SP1, CU-SP2, CU-SPR	35344	35344	CU-SP1, CU-SP2
35405	-	CU-SPR (adult/pediatric)	35405	35405	CU-SPR (adult/pediatric)
35390	-	Metrax-Primedic HeartSave Pad, AED, AED-M, 6, 6S since S.N. 739xxxxxxx	35392	35392	Metrax-Primedic HeartSave Pad, AED, AED-M, 6, 6S since S.N. 739xxxxxxx
			33388	35380	HeartSave Pad, AED, AED-M, 6, 6S up to S.N. 738xxxxxxx
			35380		Defimonitor XD1, XD3, XD10, XD30, XD100, XD300, XD110, XD330
33458	33459	Progetti Rescue Sam, Rescue 230, Rescue Life 5.7" and 7" with cable	33458	33461	Progetti Rescue 230, Rescue Life 5.7" and 7"
			33420	33428	Rescue Sam with cable
				33461	

COMPATIBLE DEFIBRILLATOR ECG THERMAL PAPER - DIFFERENT BRANDS

GIMA code	MANUFACTURER AND MODEL	Type	Width mm	Length m	Sheet	Marker	Grid*	Min. order	GIMA code	MANUFACTURER AND MODEL	Type	Width mm	Length m	Sheet	Marker	Grid*	Min. order
33099	Bexen Reanibex - 800	roll	106	20	-	no	or/red	10	32988	Physio control Life Pack 9	roll	50	30	-	no	or/red	20
33100	Corplus Corplus 3	roll	106	20	-	no	or/red	10	33032	Life Pack 8-9-10	roll	50	23	-	no	or/red	20
33103	Hewlett Packard 40457C/D,	roll	50	30	-	no	or/red	20	33033	Life Pack 11-12-15-20	roll	107	25	-	no	or/red	10
32968	Mindray Defi/Monitor	roll	50	20	-	no	or/red	20	32963	Schiller Defigard 3002	roll	50	25	-	no	or/red	20
33029	BeneHeart D3, D6	roll	50	20	-	no	or/red	20	32983	Defigard 6002	pack	80	70	300	yes	or/red	20
32988	Nihon Kohden Def. TEC 5521, 5531,	pack	50	100	200	no	or/red	25	32988	Welch Allyn Pic-30	roll	50	30	-	no	or/red	20
33102	BSM 2300, 3000	roll	50	30	-	no	or/red	20	33106	Zoll Defi/1400 (P/N 8000-0300)	pack	90	90	200	yes	or/red	20
33103	Tech-8300	roll	75	25	-	no	or/red	10	32964	X Series	roll	80	25	-	no	white	10
	Philips Defibrillator MRX	roll	50	30	-	no	or/red	20	33647	R Series, X Series	roll	80	25	-	no	or/red	10
	Heartstart XL	roll	50	30	-	no	or/red	20									

*or/red: grid could be orange or red colour

REANIMATION MASKS



• **34218 CPR MASK - pocket resuscitator**
The CPR pocket mask is a non-invasive, latex free mask with an integral non-rebreathing valve and efficient filter to prevent direct mouth contact with the patient's face when providing mouth-to-mouth resuscitation. It can be used on adult, child and infant. Supplied in a compact plastic box with a corrugated tube for an easier use.



• **34211 VENTO AIRWAYS - Poly-bag**
• **34213 VENTO AIRWAYS - Key-holder box**
Disposable mouth to mouth airway with a one way polycarbonate valve and non-toxic pvc protective foil to avoid any direct contact between patient and rescuer. Available in 2 versions: in poly-bag or key-holder box. Instructions in GB, IT, FR, DE.

DEFIBRILLATION KIT

• **34219 DEFIBRILLATOR ACCESSORIES KIT**
Pocket mask kit with accessories for reanimation and defibrillation. Kit includes:
1 pocket mask,
1 antiseptic wipes,
4 vinyl gloves,
1 disposable razor,
3 packs of sterile gauzes.



DAE-AED SIGNS


• **35326 DAE-AED ALUMINIUM FRONT SIGN 29x36 cm**
• **35327 DAE-AED DOUBLE-SIDED ALUMINIUM SIGN - 34x36 cm**
Signs in aluminium foil with thickness 5/10. Prepainted white with silkscreen in green colour. Pre-drilled for wall mounting.

35326



Signs comply with UNI EN ISO 7010:2012

35327



350 g