

SUN-III-45L Autoclave B Class



Product details

Model	Power	Voltage	Frequency	Tank Capacity	Appearance size	Inside Dimension	Outside Dimension	G"W
SUN45-III	3500W	AC220V/110V±10%	50/60Hz±1Hz	45Litres	530X770X530mm	Ø319x620mm	610×840×580mm	100kg

Technical parameters

The rated voltage	AC 220V
Frequency rate	50±1Hz
The most powerful	3500VA
Should the product	45L
Class	European Class B standard
Sterilization temperature	134°C/121°C
Special sterilization	Inactivated AIDS (HV) and hepatitis b (HBV) mad cow virus and spore
Drying procedure	Powerful vacuum drying, instrument residual humidity < 0.2%
Show	LCD display, touch the keys
Test	Bowie&Dick steam penetration test, helix spiral test, Vacuum Vacuum test
Security	The relief valve, pneumatic locking, breakdown self check system
Sterilization record	With printers, and USB interface to record the sterilization process
Water supply system	Open water tank cleaning disinfection can, with functions of water quality detection and debris filter, clean water tank volume: 5 l, wastewater tank: 3 l
Water quality detection	Whether can detect the water quality conform to the requirements of the equipment
Cleaning procedures	Automatic cleaning internal line and steam generator
The sterilization chamber configuration	Layer 3 activities tray
Inside diameter size	Φ319mm X 620mm
Net weight	113kg
Appearance of size	780X550X540mm