

seca 336

Mobile digital baby scales (III)



- Handy and light for mobile medical use.
- Comfortable cradle.
- Large, easy to read LCD display.
- Reliable results due to automatic damping system.
- Optional measuring rod seca 232.

seca 336:

The mobile baby scale with special damping system.

Modern technology can only be of use if it passes the tests of everyday life. The digital baby scale seca 336 has all that is necessary to make daily routine in the surgery and hospital so much easier in future: as this scale is not only easy to use but also mobile and durable, it meets all the prerequisites of an ideal scale for paediatrics.



Handy and lightweight for mobile medical use.

Weighing only 3.4 kg, this scale is a real lightweight in its class. The integrated side handle makes it easy to transport. Another practical feature is the flat and thus space-saving design, making it easy, for example, to store the scale away in the drawer of a baby's changing table. It can be packed away safely for transport in the optional transport case seca 428.



Easy to operate with generously sized cradle and easy to read LCD display.

The baby feels at ease and safe in the extra large, comfortable weighing tray of the seca 336. Larger babies or small children can be weighed in a sitting position. The weight is shown on the brightly illuminated LCD display and the automatic switch-off ensures powersaving operation.



The optional measuring rod permits weighing and measuring in one step.

Reliable results due to the automatic damping system.

Absolute precision in all situations: fitted with highly sensitive weighing technology and the special seca damping system, exact results are a matter of course for the seca 336. Even when babies are very restless.

Optional measuring rod seca 232.

When fitted with the optional measuring rod seca 232, the baby scale seca 336 becomes a time-saving, combined measuring and weighing station.

seca 336

Technical Data

- Capacity: 15 kg / 33 lbs
- Graduation: 10 g / 0.5 oz
- Dimensions (WxHxD): 638 x 115 x 300 mm / 25,1 x 4,5 x 11,8"
- Dimensions/tray (WxHxD): 595 x 50 x 255 mm / 23,4 x 2 x 10"
- Weight: 3.4 kg / 7.5 lbs
- Power supply: batteries / power adapter optional
- Functions: TARE, HOLD, kg/lbs switch-over, automatic switch-off,
- Optional: measuring rod seca 232, head positioner seca 419, power adapter seca 447, switch mode power adapter seca 400, carrying case seca 428
- Approval class: III