



ABPM50/PM50

Ambulatory Blood Pressure Monitor

ABPM50 and PM50 are designed to record the patient's blood pressure information for 24 hours and help the physician interpret the patient's blood pressure by analyzing the recorded data.



Technical Specifications

- Oscillometric method
- 2.4-inch color LCD display
- 24-hour continuous blood pressure measurement
- Manual blood pressure measurement
- Battery indicator
- 4800 blood pressure measurement storage capacity
- Alarm function
- Measuring ranges;
 - SYS: 40mmHg~270mmHg
 - MAP: 20mmHg~235mmHg
 - DIA: 10mmHg~215mmHg
 - HR: 30bpm~250bpm
 - SpO₂: 0%~100% (Available for PM50)
- Measurement intervals: 5, 10, 15, 20, 30, 45, 60, 90, 120 min.
- ABPM adjustment on device and viewing analysis

Physical Features

Dimension: 128 x 69 x 36 mm

Weight: <350 g (Including Batteries)

Accessories List

- | | |
|------------|---|
| Adult cuff | Carrying bag |
| Cuff hose | SpO ₂ probe (Available for PM50) |

Software Features

- Detailed analysis
- Entering demographic information
- Alarm range adjustment
- Histogram, trend table, correlation and pie chart graphs
- Report preview
- Report format adjustment
- Red marks when the limits are exceeded

Merkez Mah. Kağıthane Cad. Firuze Sok. No.5 DAP VADİ S Ofis K. 4 D. 82-83 Kağıthane - İstanbul /Turkey

+90 212 803 2808

+90 212 803 2808

info@semcares.com.tr

www.semcares.com.tr