

Dr. pen[®] A20

Dr. pen Microneedling System

FDA 510K cleared

Class II medical device



Medical Device Certification

Owner/Operator: Guangzhou Ekai Electronic Technology Co., Ltd.

Device Name: Dr. pen Microneedling System

510(k) Number: K230420

Regulatory Class: Class II

FEI Number * : 3013360665



FDA U.S. FOOD & DRUG
ADMINISTRATION

Guangzhou Ekai Electronic Technology Co., Ltd.


Re: K230420

Trade/Device Name: Dr. pen Microneedling System

Regulation Number: 21 CFR 878.4430

Regulation Name: Microneedling Device For Aesthetic Use

Regulatory Class: Class II



Advantages

LED Digital Display

1~6 gear Adjustable speed



0~2.5 mm

Adjustable needle length



Type-C USB
connector



Large capacity battery ensure
continuously long working time

700mAh



3L Connector

The connection of cartridge is with 3 locks design. The cartridge is locked the bayonet easily and firmly, the procedure will be more stable and high efficiency.



DISPOSABLE
PROTECTIVE BAG

Clean and hygienic to prevent cross infection



MEDICAL.GR


Dedicated to Health




Solve various Skin Problems

Accelerate pigment metabolism, whiten and tender skin, remove acne mark, dilute chloasma and reduce acne pits, improve stretch marks, reduce pores and fine lines.




 A20 with low battery warning to help users know the battery remaining power intuitively. Better control the working hours

 The symbol "L" will blink on LED screen for 5 times when the battery is low, and then the Dr.pen will be turned off automatically.

Intuitive low battery warning

MEDICAL.CR
Dedicated to Health



Product Parameters



Specification			
Brand	Dr pen	Model	A20
Speed Level	6 Speeds	Power Model	Wirless
Battery Capacity	700 mAh	Handle Material	Metal
Needle Adjustment Range	0 ~ 2.5mm	Available Catridges	14 Pins
510(k) Number	K230420	Adapter	US/EU/UK/AU/KR
Regulatory Class	Class II	Specification for 1 single package	20 * 16 * 8.5 cm/1 pc/ 700 g
FEI Number*	3013360665	1 carton	? cm / ? pcs / ? kg
Customization (OEM & ODM)	logo, private label, package, new shape,change color.	Packing List	<div>① Box</div> <div>② Dr pen A20</div> <div>③ 14 pins * 2 pcs</div> <div>④ Manual</div> <div>⑤ Inductive base</div> <div>⑥ Adapter</div> <div>⑦ Cable</div> <div>⑧ Protective Sleeve * 10 pcs</div>