



Meditech
Indoor AED Cabinet

Advanced technology

Indoor AED Cabinet

1 Cabinet Specifications

Type:	Indoor AED Cabinet
Model:	MDA-E13
Color:	custom
Material:	High quality steel
Use:	Store AED / first aid supplies
Size:	H40*W36*D18cm
Weight:	6KG
Logo:	custom



2 Cabinet details



Audible alarm



Security Key



Large window



Flash alert



Customized logo



Handle design with metal lock



Suitable for



Indoor AED Cabinet

Product Highlights



- | | |
|-------------------------|---|
| Alarm system | Flash alarm system |
| Alarm Specifications | Power: 9v battery
Sound: 85dB/120db
Alarm start: Once door open
Switch: ON/OFF |
| Door | Soft Handle design, Metal Lock
Rubber material for water proof |
| Has a safety lock | Configure high-quality security keys to prevent AED lost |
| Customized logo | Logo can be customized according to customer needs |
| Observation window | With plexiglass window
Observe the condition of AED instruments in the cabinet |
| Can place multiple AEDs | Applicable to Meditech AED and Most brands AED |

Indoor AED Cabinet

1 Cabinet Specifications

Type:	Indoor AED Cabinet
Model:	MDA-E12
Color:	custom
Material:	High quality steel
Use:	Store AED / first aid supplies
Size:	H40*W36*D20cm
Weight:	6KG
Logo:	custom



2 Cabinet details



Audible alarm



Security Key



Large window



Flash alert



Customized logo



Handle design with metal lock

Suitable for



Other Brands

Indoor AED Cabinet

3 Product Highlights



● Alarm system	Flash alarm system
● Alarm Specifications	Power:9v battery Sound:85dB120db Alarm start:Once door open Switch: ON/OFF
● Door	Soft Handle design ,Metal Lock Rubber material for water proof
● Has a safety lock	Configure high-quality security keys to prevent AED lost
● Customized logo	Logo can be customized according to customer needs
● Observation window	With plexiglass window Observe the condition of AED insyuments in the cabinet
● Can place multiple AEDs	Applicable to Meditech AED and Most brands AED

Indoor AED sign

AED sign-age is important as research shows it reduces the time to first shock.
Make your AED visible to the public who need it in an emergency from.



Specification •

Product name	AED Sign
Model Number	MS-3D
Material	Plastic or Aluminium or Metal
Size	H23*W17 cm
Net weight	260 g
package	60 pcs / CTN
CTN Size	40 * 36 * 26 cm
G.W.	14.8 kg
Color	Customized
Logo	Customized





Advanced Medical Technology



Meditech Equipment Co.,Ltd (Meditech Group)
Address: 140 Zhuzhou Rd, Qingdao , Shandong Province, P. R. C
Tel: (+86) 532 85832673
Fax: (+86) 532 81705331
Email: sales@meditech.cn
www.meditech.com.cn