



SKR054-ET

EMERGENCY TROLLEY



SAIKANG / Medical

Healthier
your life

SKR054-ET


FEATURES AND HIGHLIGHTS

Technical parameters

▪ External size (LxWxH)	750x475x950 mm
-------------------------	----------------

Technical configuration

▪ Defibrillator Shelf	1pc
▪ Sliding Side Shelf	1pc
▪ Sharp Container	1pc
▪ IV pole	1pc
▪ CPR Board	1pc
▪ Oxygen Bottle Holder	1pc
▪ Power Outlet	1pc
▪ Central Lock	1set
▪ File Box	1pc
▪ Trash Can	2pcs



One-piece countertop

The one-piece countertop provides more comprehensive protection; the fully guarded handle makes it more comfortable to hold me; the raised grip shape makes the operation of the side-drawing board more comfortable.

Decorative gripper

Outer beveled decorative grip, better operating comfort and more three-dimensional shape.

Heavy duty construction

Wheels assembled on the same vertical line as four columns. More stable, even the drawer open inside full of medicine/tools.

Drawer

Drawer insize big & deeper!

Power supply

The power socket is movable and detachable, easy to use, located above the CPR board to convenient the use of oxygen cylinders.

Inserted oxygen bottle holde

Inserted retractable oxygen bottle holder, belts fixed on trolley main body, space saving.

Infusion stand

IV pole inserted into main body column.

Drawer holder

Card insert from downside, easy keep clean.

Defibrillator


The plate surface is made of virgin ABS material. The position of the I,V pole and the defibrillator can be exchanged for easy use.

Sorting trash can

Green+yellow sorted trash can.



CONSTANT INNOVATION




Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical (Philippines) Company Limited

Manila, The Philippines
TEL: +63 953 044 7864

US Branch

Healthward Technology, INC

17800 CASTLETON ST STE
180.
CITY OF INDUSTRY CA 91748
TEL: +1 (603)416-3063

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
Chelyabinsk,
Universitetskaya
Naberezhnaya, Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue