



# SK-AT75077A2

## ANESTHESIA TROLLEY



# SK-AT75077A2

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Size	750 × 475 × 950 ±10 mm
▪ First and secondmall drawer, drawer front dimensions	80±5mm
▪ Interior dimensions	430×380×70±5mm
▪ Third and fourthmedium-sized drawers, drawer front dimensions	120mm ±5mm
▪ Interior dimensions	430×380×110±5mm
▪ Fifth deep drawer, drawer front dimensions	240mm ±5mm
▪ Interior dimensions	430×380×220±5mm

### Technical configuration

- ☑ The entire trolley is constructed from aluminum alloy structure and ABS material
- ☑ The ABS work surface features handles on both sides for easy maneuvering
- ☑ The work surface has a recessed design with ABS guardrails on three sides to prevent items from slipping off
- ☑ The tabletop is fitted with transparent tempered glass
- ☑ Drawer interior with 3x3 dividers, freely adjustable
- ☑ Left side:Hidden retractable auxiliary workbench, hand sanitizer holder, removable file box, storage basket
- ☑ Right side:Sharps container, ABS two-tone waste bin
- ☑ Back side:Dual-row adjustable storage boxes, 5 per row, total of 10
- ☑ Bottom: 4-inch luxury double side silent casters, 2 with brake functionality

#### Material

The entire trolley is constructed from aluminum alloy structure and ABS material.



#### Table-Board

The ABS work surface features handles on both sides for easy maneuvering;The work surface has a recessed design with ABS guardrails on three sides to prevent items from slipping off;The tabletop is fitted with transparent tempered glass.

#### Left side

The left side of the trolley is equipped with a hidden retractable sub-workbench, hand sanitizer holder, removable file box, and storage basket.



#### Right side

The right side of the trolley is equipped with a sharps container and an ABS two-tone waste bin.



#### Back

Back side: dual rows of adjustable storage boxes, 5 per row, total of 10.




#### Castor

4-inch luxury double side silent casters, 2 with brake functionality.





## **CONSTANT INNOVATION**




Design



Research



Manufacturing



Maintenance



Service

## **WORLDWIDE BRANCHES**

### South America Branch

#### **Saikang Medical Technology, Inc**

W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

#### **Saikang Medical Technology Ltd**

Mezzanine Floor, Taj Tower,  
upper hill, Nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

#### **Saikang Medical (Philippines) Company Limited**

Manila, The Philippines  
TEL: +63 953 044 7864

### US Branch

#### **Healthward Technology, INC**

17800 CASTLETON ST STE  
180.  
CITY OF INDUSTRY CA 91748  
TEL: +1 (603)416-3063

### Russia Branch

#### **LLC 'SAIKANG'**

454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya  
Naberezhnaya, Russia  
TEL: +7 (919)300-00-06

### Mid East Branch

#### **Saikang Medical Supplies (Mid East) LLC**

Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

# SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

 +86 18021231901 ·  24/7

 export@saikangmedical.com

 www.saikangmedical.com



Company video



General catalogue