



SKR-CT750-1

TREATMENT CART

Registration code: SKM-B



SAIKANG Medical

Healthier
your life

SKR-CT750-1

FEATURES AND HIGHLIGHTS

Technical parameters

▪ External size (LxWxH)	750x475x920 mm
▪ Castors	125 mm



Material

Steel and aluminum structure inside, ABS materials outside, easy to clean. Powder coated columns.




Table top

Guardrails on top of trolley is made of stainless steel, table top with soft glass, avoid possibly falling down during movement.



Push handle

With push handle on both sides, convenient to medical staff.




Drawers

Drawer with 3*3 abs inserts inside, which can be separated freely as required.



Drawer card holder

Card insert from downside, easy keep clean.



Tray

Tray in the middle.



Castor

High-quality medical silent universal castors, flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.



Trash can

Two trash cans for waste sorting.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc
W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd
Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd
Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited
DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG'
454004, Russian Federation,
Chelyabinsk,
Universitetskaya Naberezhnaya,
Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC
Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ export@saikangmedical.com
🌐 www.saikangmedical.com



Company video



General catalogue